

Product datasheet for **SC210755**

NUSAP1 (NM_001129897) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	NUSAP1
Synonyms:	ANKT; BM037; FLJ13421; LNP; NUSAP; PRO0310p1; Q0310; SAPL
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001129897
Insert Size:	911 bp
Insert Sequence:	>SC210755 3'UTR clone of NM_001129897 The sequence shown below is from the reference sequence of NM_001129897. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CGAAGGGGCTCATTTTGGCTGAAGATATAATTTTTTAACATCTTGTAATATTCCTGTATTCTCAA CTTTTTCTTTTGTAATTTTTTTTTTGTGTCATCCCACTTTAGTCACGAGATCTTTTCTGCT AACTGTTATAGTCTGTGTAGTGCCATGGGTTCTTCATGTGCTATGATCTCTGAAAAGACGTTATCAC CTTAAAGCTCAAATTCCTTGGGATGGTTTTTACTTAAGTCCATTAACAATTCAGGTTTCTAACGAGACC CATCTAAAATTCGTTTCTAGATTTTAAATGTCAAGTCCCAAGTCCCTGCTGGTTCTAATATTA ACAGAACTGCAGTCTTCTGCTAGCCAATAGCATTTACCTGATGGCAGCTAGTTATGCAAGCTTCAGGAG AATTTGAACAATAACAAGAATAGGGTAAGCTGGGATAGAAAGGCCACCTCTCACTCTCTATAGAATAT AGTAACCTTTATGAAACGGGGCCATATAGTTTGGTTATGACATCAATATTTACCTAGGTGAAATGTT TAGGCTTATGTACCTTCGTTCAAATATCCTCATGTAATTGCCATCTGTCACTCACTATATTCACAAAAA TAAAACTCTACAACTATTCTAACATTGCTTACTTAAAGCTACATAGCCCTATCGAAATGCGAGGATT AATGCTTTAATGCTTTTAGAGACAGGGTCTCACTGTGTGCCAGGCTGGTCTCAAACCTCCACCAATG TACTTCTTATTCATTTTATGGAAAAGACTAGGCTTTGCTTAGTATCATGTCCATGTTTCCTTCACCTCA GTGGAGCTTCTGAGTTTATACTGCTCAAGATCGTCATAAAATAAATTTTTTCTCATTGTATAGAAAA AAAAAAAAAAAAA ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001129897.1</u>
Summary:	NUSAP1 is a nucleolar-spindle-associated protein that plays a role in spindle microtubule organization (Raemaekers et al., 2003 [PubMed 12963707]).[supplied by OMIM, Jun 2009]
Locus ID:	51203
MW:	35