

Product datasheet for **SC210756**

NUSAP1 (NM_018454) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: NUSAP1

Synonyms: ANKT; BM037; LNP; NUSAP; PRO0310p1; Q0310; SAPL

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_018454

Insert Size: 890 bp

Insert Sequence: >SC210756 3'UTR clone of NM_018454
The sequence shown below is from the reference sequence of NM_018454. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CGAAGGGGCTCATTTTGGCTGAAGATTAATAATTTTTTAACATCTTGTAATATTCCTGTATTCTCAA
CTTTTTCTTTTGTAATTTTTTTTTTGTGTCATCCCACTTTAGTCACGAGATCTTTTTCTGCT
AACTGTTATAGTCTGTGTAGTGCCATGGGTTCTTCATGTGCTATGATCTCTGAAAAGACGTTATCAC
CTTAAAGCTCAAATTCCTTGGGATGGTTTTTACTTAAGTCCATTAACAATTCAGGTTTCTAACGAGACC
CATCTAAAATTCGTTTCTAGATTTTAAATGTCAAGTTCCCAAGTTCCCACTGCTGGTTCTAATATTA
ACAGAACTGCAGTCTTCTGCTAGCCAATAGCATTACCTGATGGCAGCTAGTTATGCAAGCTTCAGGAG
AATTTGAACAATAACAAGAATAGGGTAAGCTGGGATAGAAAGGCCACCTCTCACTCTCTATAGAATAT
AGTAACCTTTATGAAACGGGGCCATATAGTTTGGTTATGACATCAATATTTACCTAGGTGAAATGTT
TAGGCTTATGTACCTTCGTTCAAATATCCTCATGTAATTGCCATCTGTCACTCACTATATTCACAAAAA
TAAACTCTACAACCTATTCTAACATTGCTTACTTAAAGCTACATAGCCCTATCGAAATGCGAGGATT
AATGCTTTAATGCTTTTAGAGACAGGGTCTCACTGTGTGCCAGGCTGGTCTCAAACCTCCACCAATG
TACTTCTTATTCATTTTATGGAAAAGACTAGGCTTTGCTTAGTATCATGTCCATGTTTCTTCACCTCA
GTGGAGCTTCTGAGTTTATACTGCTCAAGATCGTCATAAAATAAATTTTTTCTCATTGTCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_018454.8</u>
Summary:	NUSAP1 is a nucleolar-spindle-associated protein that plays a role in spindle microtubule organization (Raemaekers et al., 2003 [PubMed 12963707]).[supplied by OMIM, Jun 2009]
Locus ID:	51203
MW:	34.1