

Product datasheet for **SC210776**

KRII (NM_023008) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: KRII

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PS100062)

ACCN: NM_023008

Insert Size: 896 bp

Insert Sequence: >SC210776 3'UTR clone of NM_023008
The sequence shown below is from the reference sequence of NM_023008. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAACAACAGGGGCCAAGAACAGCTCCTGACACAGGGAGCAGGCAGGGGCTCAGGCTCCTCTCCTC
AAATCAAGCCCTGGACAGGTCTCGCACCCACAAGTACTATCTGCTGCAGAGATCCTCTCATCTGTGGCC
AGGCACGCCTGTAATCCAGCACTTTGGGAGGCCGAGGTGAAAGGATCACTTGAGCCAGGAATTGGAG
ACCAGCGTGGGCAACTTAGCGAGACTCCATCTCTACAAAATGTGGTGACACATGCCTGTAGTCCCAGCT
ACTCAGGAGGCTGGGGAAGGAAGATCACTTAAGGCCATGAGTTTGAGGCTACAGAGAGCTATGATTGCA
CCACTGCACTCCAGCCTGGGCAACAAAGCAAGACCCTGTCTCTTGAIAAAAAAAAAAAAAAAAAATTAAAAAC
ACGTGTAGCTGTGGGTAAGGAAGTCGTCCCATTTTATGGATGTGGAACTGAGGCTGACTTAACTTCC
CGTGACTTGGAGCTCTCACCTCGAGACACCTCCAGACCCACATCCCTCAGAACGCTGCGCCAGCC
ACACCCAGATCAAGTTACACAGCTTTATTGGCCCCAGCAGCCACAGTTCCCCAAGGGGGCTTTGGG
ACAGGGAGGGCCCCGGCGTCCAGTGCAAACTACGTACGTGAGTCAGGCCCTGAGGGGGGTGCCCGGGG
TCCTCCTTCCCTACCCTGCCAGGAGGGCCCCGCTGGTGGGAGACAGGTGTATGGGCTCTGTGGGCTCTG
TCCCTGGCTGAATGGCCCCAGACTGGTTCAGCTGGGCTTGAGGGGCTGACAACAGGAAGTACCATCAT
CACTGCCAGAGCTCAAGATCCACACCCGACATGACATAAAATAAAAAAAAAATAAGTGTGTATCTTT
ACGCGTAAGCGGCCGCGCATCTAGATTTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_023008.5</u>
Summary:	This gene overlaps with the gene for cysteine endopeptidase AUT-like 4 in a head-to-tail orientation. [provided by RefSeq, Jul 2008]
Locus ID:	65095
MW:	32.6