

Product datasheet for **SC210778**

PLD5 (NM_152666) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	PLD5
Synonyms:	PLDC
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_152666
Insert Size:	1497 bp



Insert Sequence:

>SC210778 3'UTR clone of NM_152666

The sequence shown below is from the reference sequence of NM_152666. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
ACAGGCGGAAAGGATCCCCGGAACGTAAACATGATGAAGAACTGACAGGACAGCTCTGTATTTTACATA
TTTATACATAAAGGACTTGAGGAGAGAAAAACACTTTAATATGTCTCTTTTTTAGGAAAAAGCACAC
CTTATAAAAAATATTCTCTGAACAACACATAATACCTATCTAACAATATCTTAGCGTCGTGTACACTAA
AATTAAGACTTTTGTATTTGTTACCTCTGACAGAGAATGTCATTCTGGCGTAGAAGTTAGTTTTACGT
TTGCATACAAATTTAAGACTTTCATTCCGTCTCTTTCTAGTTTTAGAACTTTTCTGCTGACCTACC
CGGTCAAGTTGTTAGGAAGCTTAGCATGAGGTGGCTCTTTAACACCATTCTGAGAACTGTTGATTAGA
GAGGAAGGTGAGGATTTTGCCTTGGAGAAGAGACTGATAAAACCAAGTGTATTCTTGTAATCGTCT
CCAAGACCTACACAGTCAGCTAAGCATGAATAGCTTACATCTCATGCCAGCCTTCAACACCCTTGAAAA
ATATCCTTTTCTGTCTTGTCTATTTTCTTCTTATACCATAAGTTTTATTTGGTATCATTTTCTCTCTC
ACAAAAATGTTCTATTATATAGCATTTTCTGTATTACCTCAACATTACTTCTAACTATGATAACATCT
TTTGAAACATACAAATTTTTCAGAATGGATAAAGTGGAATAGCAGTCTTTTCACTTTTCACTTTTCA
ACTGAATCCTTTTTTCCCTACCGAATTGGTATTTAAAAAAATCACAGTCAGCAATATGTATCCTTTC
TCTCAAGCTCAGAAAAATTATGTAGCCACCTGCTAATATCCAGCAAGCTATCAAGTCTGCATTTTGGGA
ATAAAGTAGTTTTTACATTTTGTTCAGTTTATGGTATTATAAAGGAAAAAATCTCTATGCGTCAGTG
GATTTTAATTTCTTTAGAATCATTATGCTGTTAAGAGTATGCCAGCAACGTTTCTAGAGCTATTTAAT
ATACCTTCTGCTATTTAATATACCAATGCTTTTCAGAGAAATGCAATTTAAATACTGTGCTTAACATA
TTTTACTAAGAAGCAGAAATATTGAATTATTGTTCTATAACGTAATTGTGTGTTGTTTATATCTA
TACAATATATACCCTGTATTAATGGAGACCATTGCTACACTAGGAGAGGCATGGCCTCTTCAGAGCC
CACAGCAGGTGGCCGGGATGTCCAGTTCCAATGTTGTTCCCGACTCTGTGACTCAAAGCGGAGCCCTC
ACTGCTTGGCCTCGGGCTTCTCTTGGTAATGTAGACGGTAATAACTCACCTACCTCGCAAGGATGTTGT
CACGGTAGGAAATATTTATACTTGTGTTTGAATCAAAAGAACCTATATAAATGCTAAGCATTATAATA
TATAATCTTTTCTGAAATAAACTTTGTCTATCTTGATGGAATCAA
ACGCGTAAGCGGCCCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_152666.2](#)

Locus ID:

200150

MW:

58.6