

Product datasheet for **SC210779**

ZNF14 (NM_021030) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ZNF14

Synonyms: GIOT-4; KOX6

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_021030

Insert Size: 910 bp

Insert Sequence: >SC210779 3'UTR clone of NM_021030
The sequence shown below is from the reference sequence of NM_021030. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAAAGGACTCATATGGGAGAGAAAGTCTAAGAATATAAGCAACATTCTGAAGCCTTCGTTGTTCCAGC
TCCTTCTGAATACACAAAGAGTTCTACTGGACAGAAGCTGTGTGAATACAAACATTGTTCTAAAGCC
CTTAATTGTTCCAGTTTCTTTTGAGCACGAAAGAATCCCACTGGTGAGAAACAGTATGAATGTAAACA
TTGTGATTAAGCCTTCAGTTGACCTAATTCAATTTCAAAGACAGGACTCGTGCTGCAGTGAAATCTGTTT
TGTTTTTTGTTTTTTGTTTTTTTAAAAAAGGCCAAGCTTGTGGCTCATGCCTATAGTCCCAGCACTTT
GGGAGGCCAAGGTGCAAGGCAGGCAGATCATTTGAGGTCAGAAATTCGAGACCAGCGTGGTCAACATGG
TGTAACCCCATCTCTACTAAAAATTACAAAAAATGAACCAGGCATGGTAGTGCATGCCTCTAATCCCA
GCTACTCAGGAGGCTAAGGCAGGAGATCACTTGAACCTGGAGACAGAGTTGCAAGTGCAGTGCAGATCA
CGCCACTGCATTCCAGCCTGTGCGACAGAATGAGACACTGTCTCAAGAAAAAATAAAAGTGGGAAAAATC
TTCCAAGTGTCTCTAGAGTGTGAATATATGAAAGGAGTCAGAGTGGGGTGAAAGTCTATAAATGTAAG
ACATTTGGGAAAGCCTTCACACAGCCTTTTGAGCACACATGAGAATGTATACTGGAGAGAAACCCCTATA
AATATTAAGAATGTGGGAAATCTTCATCCTAGTTCTTTTGTGTTGTTGATGATACAAAAATATTTTC
ATTTAATAAAGATTGTCAGGTTATATTTAATAAGAAAGAAAGTATCATGTAACAGCCAGTAATAAAA
TTTTACAATCTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_021030.3</u>
Summary:	The protein encoded by this gene contains a zinc finger and a Kruppel-associated box (KRAB) domain. KRAB domain is known to be involved in the transcriptional repression of a number of zinc finger proteins. [provided by RefSeq, Jul 2008]
Locus ID:	7561
MW:	34.9