

Product datasheet for **SC210784**

CEACAM19 (NM_001127893) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	CEACAM19
Synonyms:	CEACM19; CEAL1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001127893
Insert Size:	895 bp
Insert Sequence:	>SC210784 3'UTR clone of NM_001127893 The sequence shown below is from the reference sequence of NM_001127893. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CCCTACTGCCAGCTGGTGCCAACTTCCATGGGTCTGGGCCAGGCCAGCCAGGGAGAAGACAAGGCC CCAGCCCTCCTCTGGGAGCCTCACACCTGAGACCAGCAGGACAAGGCCATTGGGGGCTGTGGGGCCGAT GAGGTGGACTCAGCCAAAGACTCAGCAGCACATGGGGCAGGTGCTCTGGCAGGGGGACAGGAGACTGTA ACAGGCCCAGGTCTTTGTGCAGCCCGTGAATGCACGCCCGCCTTCGGTCTGTTCCCTTCAAGCAAGCTGG CCTGGGCCATGTGCCTGTGAAAGGCAGGCTCTGGCCCCCTTCCATGCCAAAGTCCCCCAAGATCTGGAT ATCTGGGGACAAGATGGTGGCCTCAGGCCTGCCTCCAGGCAGTTGGTGGGCTCCCAACTGTCTGTCC TCAATGCCCTACCCCACTCCACTAGTGACCTCAGAGTCTTCTCCCTTAGGACAAGGCAGACACCCC ACCATGCGGGCCTCAGGTGGCAGAGAGGCCAGCCTCACAGGCCTGTGGCCCCACACACCAGTCCCAGC AAGGTGACCACGGCTGTGAGCCCTTCCCTGTTCAAGCAGGCCAGCCCTCTCAGAACCTGCTGCCA GCTGCTGGTCTTGGCCCCACCCTGAATCTTACTGAGTCCCTCTGGGCAGCAGTCCCTTCTCCACCCC ACCCAGCACCCGTCCCAATGTGGCCTCAGCTTGCTCCTCCCTTCCCCAACTATGCATTATTCAGC AATAAATGAGCCTTTGCTGTATGCCAGACCCAGTTCTAGGCTCTTGACGCCCTGGTAGAAAGCAAAACA AAAACCTCCGCTTTCCAGAGGGTCACTTGGCAGGGGAAAGAATCAGGAAATAAAATAAATAATGAA ACGCGTAAAGCGGCCGGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001127893.3</u>
Locus ID:	56971
MW:	30.8