

Product datasheet for **SC210786**

CNIH (CNIH1) (NM_005776) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	CNIH
Synonyms:	CNIH; CNIH-1; CNIL; TGAM77
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_005776
Insert Size:	886 bp
Insert Sequence:	>SC210786 3' UTR clone of NM_005776 The sequence shown below is from the reference sequence of NM_005776. The complete sequence of this clone may contain minor differences, such as SNPs. Red=Cloning site Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

GGCATGATCTATGTTTTGGTGAGCTCTTAGAACACACACAGAAGAATTGGTCCAGTTAAGTGCATGCAA
AAAGCCACCAATGAAGGGATTCTATCCAGCAAGATCCTGTCCAAGAGTAGCCTGTGGAATCTGATCAGT
TACTTTAAAAATGACTCCTTATTTTTAAATGTTCCACATTTTGGCTTGTGGAAAGACTGTTTCATA
TGTTACTCAGATAAAGATTTTAAATGGTATTACGTATAAATTAATATAAAATGATTACCTCTGGTGTT
GACAGGTTTGAAGCTTGAAGCTTCTTAAGGAACAGCCATAATCCTCTGAATGATGCATTAATTACTGACTGT
CCTAGTACATTGGAAGCTTTGTTTATAGGAAGCTTGTAGGGCTCATTTTGGTTTCATTGAAACAGTATCT
AATTATAAATAGCTGTAGATATCAGGTGCTTCTGATGAAGTGAATGTATATCTGACTAGTGGGAAAC
TTCATGGGTTTCTCATCTGTCATGTCGATGATTATATGGATACATTACAAAAATAAAAGCGGAA
TTTTCCCTTCGCTTGAATATTATCCCTGTATATTGCATGAATGAGAGATTTCCCATATTTCCATCAGAGT
AATAAATACTTGCTTTAATCTTAAGCATAAGTAACATGATATAAAATATATGCTGAATTACTTGT
GAAGAATGCATTTAAAGCTATTTAAATGTGTTTTATTGTAAGACATTACTTATTAAGAAATTGGTTA
TTATGCTTACTGTTCTAATCTGGTGGAAGGATTCTTAAGAATTGTCAGGTACTACAGATTTTCAAAA
CTGAATGAGAGAAAATTGTATAACCATCCTGCTGTTTCTTTAGTGC

ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCG

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_005776.2</u>
Summary:	Involved in the selective transport and maturation of TGF-alpha family proteins. [UniProtKB/Swiss-Prot Function]
Locus ID:	10175