

Product datasheet for **SC210795**

GDPD5 (NM_030792) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: GDPD5

Synonyms: GDE2; PP1665

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_030792

Insert Size: 897 bp

Insert Sequence: >SC210795 3'UTR clone of NM_030792
The sequence shown below is from the reference sequence of NM_030792. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AAGACCCTCATAGAGCGGAGTGGGCGTAGCTGAAGACATGTCTGTCCCACCTGTACCTGACACAGAAG
CTGGGGAGCCTAGGAGAGCTGGTGGAAAGTGTCTGAACTCGGAGTGCTCTGGGAGCGGGCTCCACAGC
CTCCTTGTGGGCTCCAGCCCCCTTGTACGCCGAGCCTCTCTTGAGGGGGACTCCCTGTCTCCTGAGGCC
CAGCTGGGCCAGGACTCCATCCTTTTCAGATGCCCTGCAGGCCTGGGGCTCCTTCTGGGAAGTATGGGG
CCTAGGGCTTGGTCCCCCTCTTCTGAGGCCCTCTCCTGTATCCCGACCTGGAAGCTTTGATGGGTCATG
GGCCATGCCATACCCCTGTGGCAATGGAGTGTGTGGATGCTCACCTGTGCCATCTGTCTCCTGTCTG
TGCCAGGAGGCACCTGAGTTCTCTGCTGTATCCTGCCCAAGGGCCTGGGCCGAGCCTCTACCTGAAG
CAACTCTGCTCTTCTGTCTCAGTCTCAAAGCACAAGGAGTTTCAGCCAGGAGGAAGCCAGCTGCAATGT
GGAGACACGTCTCTCTCCCAACCCACCTCATGCCACCGCCAACCCCTGCCCCAGGAGCGGGCCTGAG
CCACGTCTCTAGGAGCAGCTGGAGATGGCCAAAAGAGTGAGCTCAGGACTACTGGATCCCATGCCACG
GTGTCCAGCAGACCTCAAGGCAGAAGGGTCACCTAACCCAGGAGTCCACAGACTGATGTGACCTCAGGT
TCCCACATCAGTGGCCACAGGGCAGGGCCACCTGGTAGAAGTGTTCTGGATATGGCCAGGGTGGGTGT
GTGGCTAAAGTGGGCTGAACAGAGGGAACCTAGGGCCCTTGCCCAATGTGATTAAGCTGCCATCTTGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_030792.8</u>
Summary:	Glycerophosphodiester phosphodiesterases (GDPDs; EC 3.1.4.46), such as GDPD5, are involved in glycerol metabolism (Lang et al., 2008 [PubMed 17578682]).[supplied by OMIM, Jan 2010]
Locus ID:	81544
MW:	31.1