

Product datasheet for **SC210796**

ZFP64 (NM_018197) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZFP64 (NM_018197) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZFP64
Synonyms: ZNF338
ACCN: NM_018197
Insert Size: 898 bp
Insert Sequence: >SC210796 3'UTR clone of NM_018197
The sequence shown below is from the reference sequence of NM_018197. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCTGTCCC GCCGACTCCATTCCAGAT TAGTGCTTAAAAAACAAGGAGTGGGGGAAAGGAATTGAG
AAAAAGAAATCTTAAGTAGAATTCTCTAAAAGGTTTGCTCTTAATGTTTTCTTTGTTTTGTTTTGTTTT
TGAGACGGAGTCTCGCTCTGTTTCCAGGCTGGAGTGCAGTGGCGCTATCTTGGCTCACTGCAACGTCC
GCCTCCAGGTTCAAGCGATTCTCATGCCTCAGCCCTCCGAGTAGCTGGGACCACAGGTGTACGACATC
ATGACTGGCTAATTTTTGTATATTTAGTAGAGACGGGTTTCATCATGTTGAACTCCTGACCTCAAGTG
ATCTGCCACCTCAGCCTCCCAAAGTGCTGGGATTACAGGTGTGAGCCACCATGCCTGGCCGTGGTTTTG
CTCTTAATGTTTTTAAGGATGGTTGTGAATCCCCTGGCCCCATAATAAATTGTAATTTTATACTGCTT
ACTATAATTTTTTAACTGTAACAACCTTTGAGACCACCTCTGAATCGTCGATTATACTGTTGTAG
AATCTTAAATGTTACCAAGATGATTCCAATGAGGGGTTGAATTAATGCATTAAGTAGGAATTCATG
TGTTTGTTCCAACTTGATTTTCCAACCTAATAAAGGTTTCTGTCCATCTTATTACATTTGTGTAGTA
AATGGTACTTCCAGCCTCTCTTTGCCCATCTGGAATACTCCCAGAGTTTGGGGGTGTTTCATGTT
TTATACATGTAAGTCTGTTGGCATGAAGGACATTTTCTACATAATATGACATGGATACTTGACCCAAA
AAAAATGTTTAGTGCTAATGAGCAGAAAATGAATGGTCCATAATAAATTGATATCTGATTAATAATATA
ACGCGT AAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_018197.3
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	55734
MW:	33.7