

Product datasheet for **SC210798**

ZFP64 (NM_022088) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZFP64 (NM_022088) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZFP64
Synonyms: ZNF338
ACCN: NM_022088
Insert Size: 898 bp
Insert Sequence: >SC210798 3'UTR clone of NM_022088
The sequence shown below is from the reference sequence of NM_022088. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCTGTCCC GCCGACTCCATTCCAGAT TAGTGCTTAAAAAACAAGGAGTGGGGGAAAGGAATTGAG
AAAAAGAAATCTTAAGTAGAATTCTCTAAAAGTTTGGCTCTTAATGTTTTCTTTGTTTTGTTTTGTTTT
TGAGACGGAGTCTCGCTCTGTTTCCAGGCTGGAGTGCAGTGGCGCTATCTTGGCTCACTGCAACGTCC
GCCTCCAGGTTCAAGCGATTCTCATGCCTCAGCCCTCCGAGTAGCTGGGACCACAGGTGTACGACATC
ATGACTGGCTAATTTTTGTATATTTAGTAGAGACGGGTTTCATCATGTTGAACTCCTGACCTCAAGTG
ATCTGCCACCTCAGCTCCCAAAGTGCTGGGATTACAGGTGTGAGCCACCATGCCTGGCCGTGTTTTG
CTCTTAATGTTTTTAAGGATGGTTGTGAATCCCCTGGCCCCATAATAAATTGTAATTTTATACTGCTT
ACTATAATTTTTTAACACTGTAACAACTTTGAGACCACCTCTGAATCGTCGCATTATAACTGTTGTAG
AATCTTAAATGTTACCAAGATGATTCCAATGAGGGTTGGAATTAATGCATTAAGTAGGAATTCATG
TGTTTGTTCCAACTTGATTTTCCAACCTAATAAAGGTTTCTGTCCATCTTATTACATTTGTGTAGTA
AATGGTACTTCCAGCTCTCTTTGCCCATCTGGAATACTCCCAGAGTTTGGGGGTGTTTCATGTT
TTATACATGTAAGTCTGTTGGCATGAAGGACATTTTCTACATAATATGACATGGATACTTGACCCAAA
AAAAATGTTTAGTGCTAATGAGCAGAAAATGAATGGTCCATAATAAATTGATATCTGATTAATAATATA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

| | |
|--------------------|--|
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | NM_022088.5 |
| Summary: | May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] |
| Locus ID: | 55734 |
| MW: | 33.7 |