

## Product datasheet for **SC210818**

### TCP11L1 (NM\_018393) Human 3' UTR Clone

#### Product data:

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | TCP11L1 (NM_018393) Human 3' UTR Clone                                                                                                                                                                                                             |
| Symbol:                   | TCP11L1                                                                                                                                                                                                                                            |
| Synonyms:                 | dj85M6.3                                                                                                                                                                                                                                           |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_018393                                                                                                                                                                                                                                          |
| Insert Size:              | 904 bp                                                                                                                                                                                                                                             |
| Insert Sequence:          | <p>&gt;SC210818 3'UTR clone of NM_018393</p> <p>The sequence shown below is from the reference sequence of NM_018393. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATCCTGAGTAAGATCCTCGTCCGATCCAACTGTATGCACCCTACAGCAGCAGTATTACTACTAGCC
ACAGAATACCTGTTCTGTACTCTAATGTTGCATTGGAAAATGGCTATATAGTACATGTCTATTTAACAG
CACCGATTCCAAAGGGAAGAATATTGTGTATCACTGTTGAAAAGACTTGTGAGAAATCCACTGAATTC
TATTTTGAGAGATTGTATTTATGAGTGCAAGTTTACAAATCAAAGAAGCATTTTGTCTCGAGTTTTAC
AGGTAACACTCAGCTGTTGTCTCCACTCCTCCCCCATCATGTCTTCCAAGGAGGCTGGGACCTCTCTC
TTCTGCAATCTGGGTAGTTCTTTTCAGATGCCTCATGCTGAGCTGACGGCCATCCAAGCGCAAAAGAGAC
CAAGCCATGGCCTCACCTCCTGCCCTCCCTCAGACAGCCTTGTCCCTCACCCCTCCCTTCTGGTATTTT
TGTTAGAAGGGCTATATAACGCTTATTGCCATTTTCTCCCCCTATTTATTGCTTCTCCTGCCTGGC
CTAGCCAGGGCTCCCAAGCCCTTTTTCTAGAGGAGAGTTATTTGTCTCTCTCTTTTTTTTTTACT
TTCTTTTTTGTATTTAAACATTCTTTAACCTTCCCCCAAAGATACAGTTATTGGACTTTAAAAAGCTT
GTTCTTTTATACAAAATTACTCCTGCATCTGGTTGAGAGCAGTTATGACTTATGAGATCTAGTGAAAGG
AAAGCTTTGACTTTCAAGCCTTGAAAGTCAAGGAGGAAATGTATTGAAAAATAGTTTACTGTAATGTG
AACGAAGAGAAATGTCTGGAATCTTACTGAAAAATTGACCCTGAGACAATACACAGATGTGATGATTG
AATTAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

|                 |                                                                                                                                                                                                                                                                                            |
|-----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components:     | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                                             |
| RefSeq:         | <a href="#">NM_018393.4</a>                                                                                                                                                                                                                                                                |
| Locus ID:       | 55346                                                                                                                                                                                                                                                                                      |
| MW:             | 35.3                                                                                                                                                                                                                                                                                       |