

Product datasheet for **SC210819**

TMUB2 (NM_177441) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: TMUB2

Synonyms: FP2653

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_177441

Insert Size: 1095 bp

Insert Sequence: >SC210819 3'UTR clone of NM_177441
The sequence shown below is from the reference sequence of NM_177441. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTCCTAGTATTTGGGATGTATGGACGATAGGACATAGGAAGAAATGAAAGGCATGGTCTTTCTCCTT
TACGGCCTCCCCACTTTTCTGGCCAGAGCTGGGCCCAAGGGCCGGGAGGGGAGGGTGGAAAGGATGT
GATGGAAATCTCTCCATAGGACACAGGAGGCAAGTATGCGGCCTCCCCTTCTCATCCACAGGAGTACA
GATGTCCCTCCCGTGCGAGCACAACTCAGGTAGAAATGAGGATGTCTCTTCACTTTTAGGGTCC
TCTGAAGGAGTTCAAAGCTGCTGGCCAAGCTCAGTGGGAGCCTGGGCTCTGAGATTCCTCCACCTG
TGGTTCTGACTCTCCAGTGTCTGCAATGTCTGCCCCAGCACCCAGGCTGCCTGCAAGGGCAGCTC
AGCATGGCCCCAGCACAACTCCGTAGGGAGCCTGGAGTATCCTTCCATTTCTCAGCCAAATACTCATCT
TTTGAGACTGAAATCACACTGGCGGGAATGAAGATTGTGCCAGCCTTCTTTATGGGCACCTAGCCGCC
TTCACCTTCTCTCTACCCCTTAGCAGGAATAGGGTGTCTCCCTTCTTTCAAAGCACTTTGCTTGCA
TTTTATTTATTTTTTAAGAGTCCTTCATAGAGCTCAGTCAGGAAGGGGATGGGGCACCAAGCCAAGC
CCCCAGCATTGGGAGCGGCCAGGCCACAGCTGCTGCTCCCGTAGTCCTCAGGCTGTAAGCAAGAGACAG
CACTGGCCCTTGGCCAGCGTCTACCTGCCCCAACTCCAAGGACTGGGTATGGATTGCTGGGCCCTAGG
CTCTTGCTTCTGGGCTATTGGAGGGTCAGTGTCTGTGACTGAATAAAGTTCCATTTTGTGGTCTGCA
GCCTCTTTCTGTGACAGAGAATGGCTTTGCGCTGCCCCAGGAAATGGTAATGAGAGTAGGTAGGCCG
GGGCTACCAACGGGAGATGCAGTTTATTTACACCAGCAGCCATGGGGGCAGAGGGAATACACAGCGTTT
ACAAAGTTAGCTACCTGTACAGAATGGATTACATATGCAAAATAAAATCTCAAGACCA
ACGCGTAAGCGGCCCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_177441.4</u>
Locus ID:	79089
MW:	40.4