

Product datasheet for **SC210826**

C21orf58 (NM_058180) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: C21orf58

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PS100062)

ACCN: NM_058180

Insert Size: 899 bp

Insert Sequence: >SC210826 3'UTR clone of NM_058180
The sequence shown below is from the reference sequence of NM_058180. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CCCAGCCTGTGGACGCCTGGCCACCCTGAGTGTGAGTCACAGAGACCCTGGCCGGGGCACCCCTCCACC
CCCAGGCTTCTCAGGGCTGTGGCTGTGGCGGGACTATGGAAGGGAGCAGGGAGAGACCCTGCCACCA
CCCGGAGTGGCTACGCGAGTGTGGACTGCAGGCTCCTCCTGGGGAAGCTGGGCAGGCTCGCTTTCTGGT
CACGGGGCCATTCCAGGGGGCATCCCTTGTCCGGTCCCCTGCAGTGAGGGGCCTGTGAACCCACCA
GGGCAACAGCCCCCTCCAGGGACCCTCCTTTCTGTAGGGCGCGCCGGCCACCTGGGAGCCTCAGA
TCCCCCTCTCCATCACGGAGGTAAAGTTGAGGCCGTGGACGCCACCAGCCTGATGAAATAAAGATTTC
ATCCAGAAAGAGTTTAAGAAAATCTCAGCTCAGGAGCCAGAACAGAAAAGATCAATAGCTTTAGCCACA
TGGAAGTGTGAACATCCTGATGGGGAAAAATACCACAAAGTCTAATTTAATAACTAGGGACACCTGTGTC
CACAGCTGTGGCACAGACTCCCACTAAAGAAGTAGGGGTGCATGGCTCATGCCTGTAATCCAGCACTT
TAAGAGGCTGAGGCAGGAGGATCACTTGAGCTCAGGAATTCGAGACCAGCCTGGGTACATAGCAAGAC
CCCGTCTCTACAAAAATAAATTTAAAAAACTAGCCAGGTGTGGTGATATGCGCCTGTTATCTCAGCAG
CTCTGGGGGCTGAGGTGGGAGGATTGCTTGAGCCCAGGAGTTAGAAGCTGCAGTGAGCTCTGATGGAAC
CACTGCACTATAGCCTGGACCAAGAGCAGGCCCTGTCTCAAAAAATAAATAAAAAAGTGAAAGAAGC
AA
ACGCGTAAGCGGCCCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_058180.5</u>
Locus ID:	54058
MW:	32.3