

Product datasheet for **SC210849**

CCDC32 (NM_001080791) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCDC32 (NM_001080791) Human 3' UTR Clone
Symbol:	CCDC32
Synonyms:	C15orf57
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001080791
Insert Size:	902 bp
Insert Sequence:	>SC210849 3'UTR clone of NM_001080791 The sequence shown below is from the reference sequence of NM_001080791. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GGGGACAAGCCAGCAGCAGCAGAACAGTAATTACACACACACACACACACACACACACACACACACAC
ACACACACACACACACACACACACACACGCCGAGCAGCTGTCTCGGGTCCAGAGCGAGCAGCGTGGAGCTC
AGTGACAGCAGCAGGGAGAAATCCACTGAAGGAAAAACCCAAATTTCCACTCCACAAAGAAAACAGCT
GCAAGCCCCCAGGACTTACCTGGGGCTGGCATGTGTGACTGTCTCGGATGAAGTGAAGTGAAGTGAAGTGA
ACACTGGATCAAAATGCTGCTTTCTCTGTGTCTCACAGCTTGGCTGAGCTCTGTCTCTGCAGGTTAGA
AGTCTGTCTAAAGATCAAAATGTGAAAGTACTTGGAGAACTGAGGCCTCTTATGTGTAATGTGTAAGTTA
AGTGAGCCATATATTTCTTGCCTCTTCCGGACATTTCATGCTTGTGTCCCAAGCATTCCCTTGGTGAAT
TGTCACGTGAGTGGGGCCAGTAAGAGTGAAGTCTGCTCCTTGAATCCAAGCCCCATCTGGGGCTTCTCT
AACAAATCTGTAGTAAGTATACGGACTCCAGGGAGAGAGGCTGGGCTTCTTCTCTCATTTGTTCCCTTG
TGAACAAATGGGCAAAAGAAGTGTGAAATGTGGGTGTTTATGTCTGTGTATATGATTTTTTACTTC
ATGCATGGCTTCTCCTCAACTTCTCTGCACTTAAAAAGGGCCAGGTTCCAAATTAGACTTGTAAT
ATGGTGTTAGTGTGTTGACACTACTCCTGGATAGTTCCAAACATCTTCTTGTGGCAGGTTTCTGGCTG
AGCCCGAGCTTCCCTCCCTGTTTATTGTGTTTCATGATCAGTATGTGTTTCCATATAAACTTTTCTCAA
CGGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001080791.4</u>
Locus ID:	90416
MW:	33.3