

Product datasheet for **SC210850**

CCDC32 (NM_052849) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: CCDC32

Synonyms: C15orf57

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_052849

Insert Size: 902 bp

Insert Sequence: >SC210850 3'UTR clone of NM_052849
The sequence shown below is from the reference sequence of NM_052849. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GGGGACAAGCCAGCAGCAGCAGAACAGAAATTACACACACACACACACACACACACACACACACAC
ACACACACACACACACACACACACACGCCGAGCAGCTGTCTCGGGTCCAGAGCGAGCAGCGTGGAGCTC
AGTGACAGCAGCAGGGAGAAATCCACTGAAGGAAAAACCAAAATTTCCACTCCACAAAGAAAACAGCT
GCAAGCCCCCAGGGACTTACCTGGGGCTGGCATGTGTGACTGTCTCGGATGAAGTGACTGACCCAGTGC
ACACTGGATCAAAATGCTGCTTTCTCTGTGTCTCACAGCTTGGCTGAGCTCTGTCTCTGCAGGTTAGA
AGTCTGCTAAAGATCAAAATGTGAAAGTACTTGGAGAACTGAGGCCTCTTATGTGTAATGTGTAAGTTA
AGTGAGCCATATATTTCTTGCCTCTTCCGGACATTCATGCTTGTGTCCCAAGCATTCCCTTGGTGAAT
TGTCACGTGAGTGGGGCCAGTAAGAGTGAAGTCTGCTCCTTGAATCCAAGCCCCATCTGGGGCTTCTCT
AACAATCTGTAGTAAGTATACGGACTCCAGGGAGAGAGGCTGGGCTTCTTTCTCTCATTTTGTCTCTTG
TGGAACAAATGGGCAAAAGAAGTGTGAAATGTGGGTGTTTATGTCTGTGTATATGATTTTTTACTTC
ATGCATGGCTTCTCCTCAACTTCTCCTGCACCTTAAAAAGGGCCAGGTTCCAAATTAGACTTGTAAT
ATGGTGTTAGTGTGTTGACACTACTCCTGGATAGTTCCAAACATCTTCTTGTGGCAGGTTTCTTGCTG
AGCCCCGAGCTTCCCTCCCTGTTTATTGTGTTTCATGATCAGTATGTGTTTCCATATAAACTTTTCTCAA
CGGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_052849.5</u>
Locus ID:	90416
MW:	33.3