

Product datasheet for **SC210852**

FBXL20 (NM_032875) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FBXL20 (NM_032875) Human 3' UTR Clone
Symbol:	FBXL20
Synonyms:	Fbl2; Fbl20
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032875
Insert Size:	2000 bp



[View online »](#)

Insert Sequence: >SC210852 3'UTR clone of NM_032875
The sequence shown below is from the reference sequence of NM_032875. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CGCTTCTGCAGATGCTGCATCATCCTATGACAATGGAGGTGGTCAACCTTGGCGAACTGAGTATTTAAT
GACACTTCTAGAGCTACCGTGGAGTCTCTCCAGTGGAAAGCAACCCAGTGTCTGAGCAAGGGTTACAA
AGTGAGGGAGGGCAGTGTCCAGATCCCCAGAGCCACACATACATACACATACACACCTTACCCCATC
CACTCTAGCTTTGTGACCATGGGACTGAAGTTTGTGATGGCTTTTTTATCAAGTAGATTGGTAAATTT
AACCATTCCTGTTGAGGTGCCATAAGAAAATCATAGGCCAAGATAGGGAGGGGCATTCCAGCAAACCC
CGTGTTAATGCTACTGTGGTTTTAAATTTTTGTCTAGGGGTTTCTTTGGGGATTTAGAACAGCATCT
GCTGTCCTCCGGGTCAAGAAAAGCATGGAAAGACAATATATGATGTACCCAGGGACCAGAAAGAAAAT
TTCTTTGCATCTTAGAAATGGTAGACATTCATTGTGACTAAAGAGCTTCTATGCTTCCTTGTTCATG
CCAACATGCTGAGCATGCTCACAAAGAAGGCTCGTCCATTCTCTCTGTGTTTTAGATTTGGCCAGAG
GTTTCCTAAATGGTTGCCTTGAAATCACTGTGGTCCAAATGTAATTCTTACACACTCAAATTATCACTG
TCTGTAGCACACTTGTGCACCTGTCTTACATTCTCTGTTGCTCCCCCCCACACTCTTGCTCAGTCTGTC
ACCTGTTCACTGTGCTTACTCACTCAATTGTTACCTTTTTGCTGTTGCTGTTTACAGTTTGCATTTT
GAATGATTAGTTGGGATTACCAACATTTTTTAAAAAGATATTATCAATAAATATTTTTTAATTCTAA
ATTTTACCATTAGGGGACCCTGCCTGAATCTTTGCCAACCTGCAAGGAACTCTAATGAAAGCAAGAAA
AGTTATGAGAAGAACTGTTTTGATGAAAACAGATGCTACCTTTTGGTTTCTATTTTATCTTCTGATAT
GTGATAATTTCGCATGGTATCCCCAATATTCCCCCTACATATTCTCTTCTGTCAGACTCCACCAACTG
TATCATAAAGGAAAGAGAGAGATGTATTATTAGGAAAGAGTGAAAATGATAAATCACAAAACAGAGCCTG
AAAGTCCCAACACCCAAGTTCAGACATTAATGATTAATGTTATTATATTTTACTGAAATCTAGGAGT
GATCAGTGGCATCTCTGACCTATGAGGTGTTTTAACAGGAATGCCACTAGAGGAGCTGTGATTCTAG
CCAGTCTGTTTCTCTCTTTTCCCCCATTTACGCTTTTACGGGATACAAGAATGGCTCTATAACAGCTGG
TCAATTGTATGTCAGTGTAAATGGATTCAATCTAGTGTGGCACCCCAAGGGCCAGAACCTTCCCAT
AACCTGTCAGTCATCTCTCTTTTCTGCCCCAGTCTTTTTTTTTTTTTTTTTTTTTTTTTTAAAGACGGGT
CTCACTCTTGCCAGACTGGAGTGCAGTGGTGCATCTCGGCTCACCGCAACCTCCACTTCCAGGCTC
AAGCGATTCTCTTGCCTCAGCCTACCGAGTAGCTGGGATTACAGGCATGCACCACTACTGCCTGGCTAA
TTTTTGTATTTTGTAGTAGAGATGGGGTTTACATGTTGGCCAGCCTGGTCTTGAACCTCCTGACCTCAA
TGATCCACCTGCCTCAGCTCCCAAGTCTGGGACTACAGGCATGAGCCACCATGCCCGGCCCTGTTC
AGACTTCTGAATGAAGTATGAATGTAGCCAGGTGAAGCATTGAGGCAGCTGAAAGAGAACTGAGAGG
AAGTAAGGGGTTAGTGTCTCTGCTTAGCCTGGGCACTAGAACTGAGGGAGGAATTTAGTGTCCAGGAC
TCAAGAGTCCAAGAAGAAGCACCTCTGTCTTCCCAAAAAATCAGAGGATAAAGGAAGTAAATGAGT
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_032875.3](#)

Summary:	Members of the F-box protein family, such as FBXL20, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]
Locus ID:	84961
MW:	77.2