

Product datasheet for **SC210854**

SLC43A3 (NM_017611) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	SLC43A3
Synonyms:	EEG1; ENBT1; FOAP-13; PRO1659; SEEEG-1
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_017611
Insert Size:	868 bp
Insert Sequence:	>SC210854 3'UTR clone of NM_017611 The sequence shown below is from the reference sequence of NM_017611. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC TGGAAAGAAAGTCCCTCTGCAATTGCAAGTTCAGAAGCCCTCACTTTTCAGCCCCGAGGATGGTTTTG TTCATCTTCCACCACCTTTGAGGACCTCGTGTCCCAAAGACTTTGCCTATCCAGCAAAACACACACA CACACACACACACACAAAAATAAGACACACAAGGACGTCTGCGCAGCAAGAAAGAATCTCAGTTGC CAAGCAGATTGATATCACACAGACTCAAAGCAAAGGCATGTGGAACCTCTTTATTTCAAACAGAAGTG TCTCCTTGCACTTAGCCTTGGCAGACCTTGACTCCAGGGGAGATGACCTGGGGGAGGAAGTGTGTCAA CTATTTCTTTAGGCCTGTTTGGCTCCGAAGCCTATATGTGCCTGGATCCTCTGCCACGGGTTAAATTTT CAGGTGAAGAGTGAGGTTGTATGGCCTCAGCTATGCTTCCTGGCTCTCCCTCAAGAGTGACGCTTGG CTAGAGAACTCACAGCTCTGGGAAAAAGAGGAGCAGACAGGGTTCCCTGGGCCAGTCTCAGCCCAGCC ACTGATGCTGGATGACCTTGGCCTGACCCTGGTCTGGTCTCAGAATCACTTTTCCCATCTGTAAATTTG AGATGAATTTTGGTGTGAAAGTTCTTCTGGAGCAGATGCTCTAGAAGGTTTTAGGAATAGTGACAGA GTCAGGCCACCCCAAGGGCCATGGGAGCCAGCTGACCTGCTTGACCGAAGGATTTCTGACAGACTATCT TTGGGGATGTTTTCAAGAAGGGATATAAGTTATTTACTTTGGGCATTTAAAGAAAAATTTCTCTCGGGA ATAATTTTATAGAAAAATAAGCTTCTGTGTCTAAGGCAA ACGCGTAAGCGGCCCGGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_017611.3</u>
Summary:	Putative transporter.[UniProtKB/Swiss-Prot Function]
Locus ID:	29015
MW:	32.7