

Product datasheet for **SC210855**

CNRIP1 (NM_015463) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	CNRIP1
Synonyms:	C2orf32; CRIP-1; CRIP1
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_015463
Insert Size:	849 bp

Insert Sequence: >SC210855 3' UTR clone of NM_015463
The sequence shown below is from the reference sequence of NM_015463. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site **Blue**=Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

GGTGAACAAGGAGTCTTCCTCTGAAGTGGCTCTTTCTGATGCTACTGGACAATCTTTTCAGAAATCTA
CTTTTAGATAAACCAGCGCAGGCCTTAAACAAGGCATGCCACACCATTTGCTGTTCCCATCTGGTACTAA
TGACCTACTAGCCCTGGTTATTTTGAAAGTGTAATTCCTCTGACACAATATCCCTCACATACAAGATG
GTGATGTACCACTCCCGAGAATATTCCTCTGTATTTGTGCTGATGTTCTGCATCCAGTTGTTCTCG
GTCATACTTTTGACCACTTGTGACTGGAGTTCAGTGGCCCTGGCAGGCTTGCTCCTGCTCTTGACCATTC
ACTGACTAACTTTGGTGTGTTTGTTCCTCAAGTTAAGTGATTCCTCTTTTGTTCATGTTAAATTTA
AAAATAACAATGTGTATGGGTCCTCCCATGTGTAATATGGTAACATGTAACCTGCAGTGTGGCAGCTT
TCAAAGCAGGCTTTGTGAAAAATGTAATACAAACAGCAGTGAATGGGACTCAAATGTTGTGCTTCCTATAA
ACAGCTCCGCTCTTTCAGGGAAGGATGGTAACAACTAGAAGGACAAATATGTACGTATTTATAACGTAT
TAAACTCTTTAAGTAGCTTAAGGTATTGTGCAATGGCCTAGCCTAGTAGAAATGGGGAAAAGCATTG
CTGTGGACCATTTGTTAAAGTGACAGGAGTTGTAGGGTTACCCCTTTGACAAGCTTCCATAGCTTCAGAC
ACGCACATTGATGCATCCCTACAGTCTTCTCACTTGAACCTTTGTAACCTGCATATGGACCTGTAATC
AACTTTTGA

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCG

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_015463.2</u>
Summary:	This gene encodes a protein that interacts with the C-terminal tail of cannabinoid receptor 1. Two transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2013]
Locus ID:	25927