

Product datasheet for **SC210874**

Collagen VI (COL6A2) (NM_058175) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Collagen VI (COL6A2) (NM_058175) Human 3' UTR Clone
Symbol:	Collagen VI
Synonyms:	BTHLM1; PP3610; UCMD1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_058175
Insert Size:	903 bp
Insert Sequence:	>SC210874 3'UTR clone of NM_058175 The sequence shown below is from the reference sequence of NM_058175. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACAGGTTTGGACGGAGCTGTTTTGTGCTTGAAAGGTTTTCTCGGGTCCGTGGTGTCCCCAAAGGTGCC
ACCGTGCGGGTCTCCTAGCTCCCTGCCAGCTTCTGTCCCTGTGCTCACTGCCCCACGCCTCCTGCCA
AGGCCGAGCCACACACCCGCTCCACCTGCATTTCTCTACCGACTCGCCAGCCAAATGCCGCTTCTCA
CTCTGGCCTCGCTGAGCGGCTGCCCGAGGAGAGCTCTAGGCCGACGCCACCGCAGGCCTTACAGTCT
TCTCTGGACGCTCCCTTGACATGCACCGTGGCCTGGCGGCGAGCCCCGGTCACCTTCTCCGCACGG
AAGAGGGGCCGGACGCCACCTTCCCCAGGACCATTCCCCTGATCCAACAGTTGCTAAACGCCACGGAGC
TCACGCAGGACCCGGCCGCTACTCCCAGCTGGTGGCCGTGCTGGTCTACACCGCGAGCGGGCCAAGT
TCGCCACCGGGGTAGAGCGGCAGGACTGGATGGAGCTGTTCAATTGACACCTTTAAGCTGGTGACAGGG
ACATCGTGGGGGACCCGAGACCGCGCTGGCCCTCTGCTAAAGCCGGGCACCCGCCAGCCGGGCTGG
GCCCTCCCTGCCACACTAGCTTCCCAGGGCTGCCCGGACAGGCTGGCTCTCAGTGGAGGCCAGAGATC
TGGTGGCCACCGTGTCCCTTGCTGCGGCTGCATTTCCAGTCTCTCCTCCGTCTTCTGTGGCCGCTCT
CTTTATAAGAACCCTGGTCATTGAATTTAAGGCCACCCCAAGTCCAGAATGACCTCGCAAGACCCTTA
ACTCAC
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_058175.3</u>
Summary:	This gene encodes one of the three alpha chains of type VI collagen, a beaded filament collagen found in most connective tissues. The product of this gene contains several domains similar to von Willebrand Factor type A domains. These domains have been shown to bind extracellular matrix proteins, an interaction that explains the importance of this collagen in organizing matrix components. Mutations in this gene are associated with Bethlem myopathy and Ullrich scleroatonic muscular dystrophy. Three transcript variants have been identified for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	1292
MW:	32.1