

Product datasheet for **SC210883**

GPR17 (NM_001161416) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: GPR17 (NM_001161416) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: GPR17
ACCN: NM_001161416
Insert Size: 914 bp
Insert Sequence: >SC210883 3'UTR clone of NM_001161416

The sequence shown below is from the reference sequence of NM_001161416. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCTCGCTGAGTGCCAAGTCAGAGCTGTAGCGGGGGCGCCGTCCAGGCCGAGCGCAGACTGTTTAGG
ACTCAGCAGACCCAGCAAGAGGCATCTGCCCTTTCCCGAGCCACCTCCCGAGCAAGCAACCTGAAATCT
CAGCAGATGCCACCATTCTCTAGATCGCCTAGTCTCAACCCATAAAAAGGAAGAACTGACAAAGGGG
ATCCATCGGCCACCCCTCTGCAGGGCTTGTGATGGCTACAATGGCTCCTAGACACTCAACGACTTCAT
CTGTGGCAGGGAGAGAGGAGGCCGGAAGAACAACCCCTGAACAATGGAGGCCCTTTCTTTCCCGCTAGGC
TCCCGAGCTCCTTTCCCGCTACAGAATCGCTCATCGGCGAGGCTCAGCAGAAAGACCTGAAGGCAGGCT
GCAAAATGACCCAGAAGAGGGACCTGGGAGTCTGGTGGGGACGGGGAGGGAGTCTCAATACTCCTTTGC
AGTGCAAGGTAAGTCTGAGTCCCCTCTGTAGTGCCTTCCAGACACACACTGCCTGAGTTGAAGAGACA
CAGGCCACACATTTAGGCTGGTTGCCAGCGGACGTGAGCACTCACGGCTGCAGGGACTCAGCACAGC
TCTGGATTCTGGATCTCTCTGCTGTAACCCACGCACAAGCCTGCAACCCCGAGAGCTCTTTGACAGG
CTCCAGGCCTCCAGTCTGGACAAGCATGTGCAGTACGGGAGCTCAGCTCAGGCCAGGGCTGGGCT
GTGCACCTGCCTCCACTGACCCAGACCACTTCTCCAGAGAGGCTCTCTCCGCTGAGCTATTTCC
CTTGCTAGTGTGAGATATTTCCCTAACATGCTCTTTTTTGTATTTGTTGTACGGACCATAAATATAA
CTGTAGCTTTAAGACTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001161416.2
Summary:	Dual specificity receptor for uracil nucleotides and cysteinyl leukotrienes (CysLTs). Signals through G(i) and inhibition of adenylyl cyclase. May mediate brain damage by nucleotides and CysLTs following ischemia.[UniProtKB/Swiss-Prot Function]
Locus ID:	2840
MW:	33.3