

Product datasheet for **SC210888**

ZNF536 (NM_014717) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF536 (NM_014717) Human 3' UTR Clone
Symbol:	ZNF536
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014717
Insert Size:	925 bp
Insert Sequence:	>SC210888 3'UTR clone of NM_014717 The sequence shown below is from the reference sequence of NM_014717. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
ATCAAGAGGCCAGACTTGTGTGGTAAGTGACCACTCCCTGTCCTAGTCGGTCTATCTGGACTTGCCCTTG
TCTGTTTCGTGGTCCGTCGGTGGTTATCTGCAGCTTGTTAATCGTGTAAGTCAAGAGAAGAATGTATACA
CATATGTGTGTTGAATAATTACTATTGGCATAGGTATGTGTATACACACGGTGCACCAATCTACAGTAT
ATATAGCAGAGAATCAGAGGCTAAAAATATTACCCCATATGTTCCAGTATTAGTCATGGATTGCAAAGG
CTTAGTAACTTGAGCAGGAGAGAGAAAACTCCCTCAAAGTCATAAATCCTGAGTGACAACCTGCTGGAT
GACAGATCCCTTCACCTGTGGACAACCTGGCTGGGGGTGGGGGGCTGTTCCACCAGCTCACCTGAGCAT
GTAGAGGTGGGTCTGCAGTGGTCTCGTGGGTATTACTGCTTGTGTCTGATTGTCCTGTATTTTGTAAAC
ACTTTAGAAGAATACAGAAAAGTGCAGTAATTCTTTCTCCATAGTATTTAAGCAGAAATATTGCTAG
TTTAATATTGTGTCAGGTCGTCTATTAAACCAGGAGCAGATGACAGTAAAATTTCAAGTGAATAGCACCT
TGACATCTACAACTTAAAAATGGTGATTGAAGCAAAATATGTAACTTGTACGGGGTGATCGTGTGCTT
TGGAACAGAGTATTGTTGAAGTAATTAGAAGATATATTAAGGTGTTCTGGTAATGAAGGCATGTAAGT
TATAATAATTGTAGCTTTCTGAATAAGTGTCAAATATATCTTTAAGTGTGCTGTATGCTGAGTTACAA
GTTAGGTCATTTATGAATGGAATGTAAATAATACTAAAAATGCTTCAATAACTTATCTTGGTATTGCT
AATAAAAAAAAAAAGCTGTGAAACATTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_014717.3</u>
Summary:	The protein encoded by this gene is a highly conserved zinc finger protein. The encoded protein is most abundant in brain, where it negatively regulates neuronal differentiation. [provided by RefSeq, Sep 2015]
Locus ID:	9745
MW:	35.4