

Product datasheet for SC210891

OAF (NM_178507) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	OAF (NM_178507) Human 3' UTR Clone
Symbol:	OAF
Synonyms:	NS5ATP13TP2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_178507
Insert Size:	1222 bp
Insert Sequence:	>SC210891 3'UTR clone of NM_178507 The sequence shown below is from the reference sequence of NM_178507. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CTCTGGGATGAGGATCCCTACCCAGGCAGGTGGGAGCAACCTGGCGGGTGGCTGCTCTGGGCCCACT
GCTCTTACCAGCCACTAGAGGGGGTGGCAACCCACCTGAGGCCTTATTTCCCTCCCTCCCACTCC
CCTGGCCCTAGAGCCTGGGCCCTCTGGCCCATCTCACATGACTGTGAAGGGGTGTGGCATGGCAGG
GGGTCTCATGAAGGCACCCCATTTCCACCTGTGCCTTCCTTGCGGGCAGAGAGGGAGAGAAGGGCTC
CCCAGATCTACACCCCTCCCTCCTGCATCTCCCTGGAGTGTTCACTTGCAAGCTGCCAAAACATGATG
GCCTCTGGTTGTTCTGTTGAACCTCTTGAACGTTTAGACCTAAAAGGAGTCTATACCTGGACACCCAC
CTCCCCAGACACAACCTCCCTTCCCCATGCACACATCTGGAAGGAGCTGGCCCTCAGTCCCTTCTACT
CCCCAACAAGGGGCTCACTATCCCCAAGAAGGAGCTGTTGGGGACCCACGACGAGCCCTGTACTGG
ATTACAGCATATTCTCATCTCTGGCCCCGAGGCTGCCTGTGGGGCGAGTGAGACCTCCCATCACTGAG
ACAGATCACAGACCAGAGTGCCCTTCCCGACCTGGACGTTGCCTCCAGAGCAGGCACCACTCTTTC
CCTCTCTACACAGAAATATTTTGTAAAGTTCTGGGGCAGGGAGGGAGCATGAAGTACGAGGAAAACCT
GAATCCAGATTTTGTGCAAGTATTTATCATTCTACCAGAAATAAACGTTTTAAGTTTTACTTG
ACTAATGAGACCCAGAGTTTGGAGAAAACCTTTGGCCAATGCTGCCACCTGATGTCAGAAAAGTGTCCTC
ACACCCTAGCAGTGCCCTATCTTGAACAAGAACCTTCAAAGCACCTACTGTGTGCTCAGCCATTTGAG
GAAGGAAGGAGGAGAAGGAAGATGTTACTAGGGAAGGATGAGATAAACTTCTGCACCCAAGACAATGA
GACAGACATAGCTGCAACCGTAGTAAGCCAGTCAGAAATAGCCAGCGCAAGGCAAGAGATGGGGTGGA
GATTGGAACCCCGCTTCAGATCTGGGCTCGGCTACTTACCTGCTGTGCAGCCATGGGTCAAGTTGCTTG
ACCTCTGTGTGCCTCCCTTAGCTATAAAATGAGCTTACTTAAGA
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_178507.4</u>
Locus ID:	220323
MW:	45.4