

Product datasheet for SC210897

TRRAP (NM_003496) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TRRAP (NM_003496) Human 3' UTR Clone
Symbol:	TRRAP
Synonyms:	DEDDFA; DFNA75; PAF350/400; PAF400; STAF40; TR-AP; Tra1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_003496
Insert Size:	916 bp
Insert Sequence:	>SC210897 3'UTR clone of NM_003496 The sequence shown below is from the reference sequence of NM_003496. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATGGACCCCGCTGGCACCCCTGGCTGTGACTGTGGCCGCCACGGCCACCCGGAATGTGAAGGGCGCTC
CGGGCTCTGAGCCCGCAGCTTTTACGACTTCTCCCTGCCTCGTTCCTTATATTCACAGAAGCCCCATAG
TTTCACTGGGTTGCGTTATTTTCTGGTAGTTTGCCTGTAAGAAAGGGAGAATATAGTTTTAGAGGAA
GCTGAAGTATGACGATGCTGGGCGAAGCGGTTGGAAATGGCAGAGCTGAACTTATTCCAAGCTTTCAA
AATAATCTTTTAAAGAAGCCAGGATTCTCCGGTCTGGAATTTCTGAGTGAGTCCTTTTTTTATGGTGCC
TCCCTCTGTGAATGTACAGGCGGAAGTGTACGAACAGCTCCCTTCCATCCATTTTTAACTCTTTCGGAA
ATAACACCTCACAGCAGCTTCGTGCTTTTGTACAGACCTTTGTAACAAGTGTACAGAAAACCTATTTTG
TTTGAGAAACAGGAGTTGATGAACCCATCATGCTGGTTTTTCTCTGAGCAGAAAGTTTAGGCTGTACA
CAGCCAGCCTTGGGAATCTCGTTGAGCGTTGCGGTGGATCCACGGGGCCAGGCCACCTGCGGGAGCG
CCACACGCATCCACTTCGGATTCAGTGGGTGAAGACAGAACTCTGAGAGTCTGCAGGCGGCTCCTGTGC
TTTTTATTTCTGGCTCTTCGGATGTCTTCTAGACATTTACTATCACTGCACCTGAAGAAAAAATCACTT
TTACCTTCCTAATTTAAAAAGACAAAACAGAAATGTACGTTTCCTTCGCTAGCTTTAGTCTTTCTGTTCC
CATTTTTATAAATCTGAGCATTGATAATGTTCTATCTAAATTTGTACAGTGTGATTTTTTTTTTAGAA
TAAATATTTTATAAAAGG
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_003496.4
Summary:	This gene encodes a large multidomain protein of the phosphoinositide 3-kinase-related kinases (PIKK) family. The encoded protein is a common component of many histone acetyltransferase (HAT) complexes and plays a role in transcription and DNA repair by recruiting HAT complexes to chromatin. Deregulation of this gene may play a role in several types of cancer including glioblastoma multiforme. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2011]
Locus ID:	8295
MW:	35.2