

Product datasheet for **SC210903**

KCTD14 (NM_023930) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	KCTD14 (NM_023930) Human 3' UTR Clone
Symbol:	KCTD14
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_023930
Insert Size:	908 bp
Insert Sequence:	>SC210903 3'UTR clone of NM_023930 The sequence shown below is from the reference sequence of NM_023930. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
ATTTATTCATTCACCTTCACCTGGTGGTGCATCCTCAGGAGCAGAGACTGTTATGAATTCTGGCGTGGCT
TATGAAATTTAAAGTTGCCATCAAAGCCATTTCTTTTAAATTTTCAAAACATCAGGCAATTTCCAGGT
TGGTCTAGAGTCTTGCCACTAAATATTGATCACTCGTTTAAAGACTTTCCACTCCATTGCAACTGATGC
CACTATATTTGCCTAGCAACTTGCACTACTTCTTTTCAAAGCCTCATGTATCTCCCAGACCCTTCTC
TTGAAGTCCAATAACAAGACCAAGTAAGAATGTTTCAACAATGCGTTGGCAAGAGATGTGAGATGACAA
CAGGAACATACAAGATACTGTGAATCTAGATGTTCTGACCTAAAGATGTAGTCTACATAGCCCCAGCTT
GGGGTCCAATCCATCTGTCCCTGGCATGTGCCTTCATGTAGTAGGTGCTTTCCTGATCCCCTTTGCGAG
ATGCTGTGGGTGCTAACACCTCAGAGCTGTCTCTCTAGAGTGGAGGTTTTCAAAGTGCATCATCA
GCATTACCTGTGAAGTGTGCTGGAATAACAATCCTCAGGCCCCACCTCAGACCTACTGAATCAGAATCT
CTGGGGTTGGCACAGCATTCTGATTTACCAAACCCTCCAAGTATTTGATGTATTCTAATTTTGA
CCATCTCTAGAAAAGAATTGCTACCTCTTGATGGAGGTACAAAAGACTGACCTTTACATCAAGGAAC
TTCCTTTCCAGAGCTCCTCATGGAATCAAGCTGAAGTCAAGTCTTCTTCTGAGAGCACATTTACTCA
GTTTTTTTCTCTGTCTACGCTGCTTCCCTCACTCCCCTTCTCCTAAGAGCACTCCATCAATAAACCA
CTTGACAGAGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: Sgfl-Mlul



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_023930.4</u>
Locus ID:	65987
MW:	34.1