

Product datasheet for **SC210925**

CTTNBP2 (NM_033427) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	CTTNBP2
Synonyms:	C7orf8; CORTBP2; Orf4
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_033427
Insert Size:	916 bp

Insert Sequence: >SC210925 3'UTR clone of NM_033427
The sequence shown below is from the reference sequence of NM_033427. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AATGAACACCTAGAAAAACCTAACAAATAGGCCTGCCTACAATATTCTCATTATTAACCTCCTCTATT
CACACAGAAACCAAGGACAACAAGATACAAATACCTTTAAATTGATAAAATGTTTTCACTGTAATATAA
ATTATGAGTGCGGGACCACTATTCAATTTTTTTTTGTAAAGAATGATTTATAACCAACATTTTGTTC
TTTTTTTTTTAATTCTGTAGGATTGAACCTCAATACCAGGAAGTTATGCCAACTTTTCAAGTAGTTTTA
TAATAGAAAAATTGAATAGTTTTAATCAAGATGTTTCAGGGGTATTAACAGAACTGTGTACTGTATA
TCTGTATATATGTATATATGTATATTTAAGGAAGTCCAGAAATCATACTGCATACATGCTACAAA
AAAGTAGCAAAAAAGTACAAAGTTGGCTTAGATAACTTGTACCGTTTTAAATGGAGATCATAATTTTAT
GTTCACCTTGATGGGAATGACTAGAATGGGCAATAACAAATCATTCTAAACAGTTATTTTCCATTAA
AAGCCATATTAAGTATTTTGTCTGTTAAATTCTAGTTTTATGTTTTAAATCTCTAGTGGTTACATTGC
CATTTGGTAGCATCATGTATATTTATGCTTATAAGTTTGCATTATGACTCCTGTAATATGAAATT
ATAATGCAATAAAAGTGTATCTTTCTTACAAAAGTTTGTACCTAAGTGGATATAGATAACTGCATAA
ATGAATAAATTCAGAAATTAACCTCTTCAATATTTCTAGTAGCACAACTGATAATATTGATATTTTCTAC
AAATGTTACTTCTAAAATTACAGACAAAATTACTTACACAATTCAGTTAAGATGCTTTGACACTGAATA
AAATATTGATTTTCATGAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_033427.3</u>
Summary:	This gene encodes a protein with six ankyrin repeats and several proline-rich regions. A similar gene in rat interacts with a central regulator of the actin cytoskeleton. [provided by RefSeq, Jul 2008]
Locus ID:	83992
MW:	36.5