

Product datasheet for SC210931

CCDC82 (NM_024725) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCDC82 (NM_024725) Human 3' UTR Clone
Symbol:	CCDC82
Synonyms:	HSPC048
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024725
Insert Size:	903 bp
Insert Sequence:	>SC210931 3'UTR clone of NM_024725 The sequence shown below is from the reference sequence of NM_024725. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TACTTTCAAGAAGAGAAGTTTGAGTTGACACATTTCATCAGAGAAGATTTTTAAATTCCTGTAAA
TGTGAAGATCATGATTCTTGTCTGTATCATGTGACATGTTGTACATTTTATCTATATCTTCATG
GCAAAATATTTGTTTAAACATGATTATCTTAATTCCTGTGAAGTGGCACTCTACACTCTAATTAAC
TTTTATCTGATGTACATAAAAGCATTATCTTAATTTTAAAGTCTGTAAATATATTTTAAAGTTATATAA
TAATGGCTTTATAACAGATGACTGTCAAGTGAATGAGCTGTTGATATCCTGTCAGTTAGTCAAAATAT
ATTGTATCTTAAATGTATTTAGACTAACGTCTACATGTATTTTACGGAATTCCTGTCCAGATCTGT
TTATTCCTTTACAGTATTAATGATTTTCACTAATATATTTTACTGCTATCATCTAAATCAGTGGGC
CATGGAAGAAACCTCATATATCATTCAAGTTCAATCCCATTGCTTCAGAATCAGATCATAGCATATCTAT
TCTAAGTAGAGGAATGCTTTTCTTTCTTGAAAGTTTCTCTTCAAATAGATTTAGGCCTCTAGA
GTCTTTCAACCCTCACATTCAGGAATAATTTTATGTAAATTTTATGCTTTTATAATGTTCTTACTTTTT
TCTATTCAATTTTGTCTATATATTGATGATTAAGATGTATTTATTTATTTTATCCATAGCTTTTTCC
ATATAAGTATGTATCTTAGGGTCAGAACTGCTGAAAGAGACAACTCAGCCAAAACACTTGGAAGCA
TATTTTGGTATCTGCATTGCTTTGCATATCTAAATTTCCATAAAGTAAATAAAGTAAATGGTTGCAG
AGTGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_024725.4</u>
Locus ID:	79780
MW:	35.6