

Product datasheet for **SC210940**

USP35 (NM_020798) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	USP35 (NM_020798) Human 3' UTR Clone
Symbol:	USP35
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020798
Insert Size:	1449 bp



[View online »](#)

Insert Sequence:	>SC210940 3'UTR clone of NM_020798 The sequence shown below is from the reference sequence of NM_020798. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGCGATCGCC GGTGGTGACTTCCACAGACTGGTCTTCTAATGTGAACCTGCTGCCAACCTGACCCCTTCCCTCCAGGAG CCAGGTAGGGCTGAGGGAAGCTGTGGAGGCAGGCCCTACCAAGAGGAAGGATGGTACAGCTCATGGCA CCTTAGTCCTCAGCCTGATGAAGGGTACACAGAGATTCTCTCAGATATGGAAGTAAGACCTAAGTCCCT TTCATTGGGGATCAGTCCATTAAACTTTACACCAAGTGCTGTTAACTTGAAGCAGCCGAGATG GGCACACACGGGTCTTTGCCTCCCCCTCTTCCCTAGCAGGCTCCCCATGCGGGAAGATCTGATGATGT TCAGGAAACAGGCTAGACCTCAGCTCCAATGTTTTGACATCAAGTACTATTTTCTTCCGACTGCTGTA CGGTATAAAGCACAGCAGGATCCAAGCCTTGCACAAAGGGGTGGGGGGGCGAGTGTCTCCTCTGGCTGT CCTGTTTTGTTTCTCATATGGGGGTGGGGGTACCTGCACTGTCTGTACCTTTCTGAGAAGAACAG AGACCGAGACCTGCCCCCTTACCAAGCGCCACTGCATGGTTTTGGGGGGGGGGCGGGGGCTAGCTTC TCACAGCAGGAGGCCCTTGCCTCCACAGCCCTCCACGCCTGCCTCAGGGCCTGAGAAGCCACCACTTG TATCCCCCTTGTGGTTAGAGTCTGATTTTACTGCAAAGGTGTTTCATGTTCTTGTGAAGTGTGGGCTC TTAGGAAGCCTGTGGGCTCCTCTGAGCAGTTGGCCTTTGTAGTGCAACAGCAGCCACCTGCAGGTTGG GTGAAGTGCCCTGACACTGCTGTAGCCCCCTTCTAAGTCTAACCAGAGCAAGACAGACACCCATGTT CATAAATAAATAAAGTAAGCCTAAGCAAGAGATTGTTCTTCTTGGACTCAGGGGCTGTGATGGCCAC TGGGTTTTGCTCACGGGTCTGGGGAGGCCAGGACTCAGCCTGCCCCCAGTGACGTGTGAGATGCCAC AGGGAGCCAAGAGCCCTCTACTGAATCTACTGAAGTGGGTGTAAGTAACGTGTGGACTGACTGACTT ACTTACCTTACTGAGGGCTGGGTGATGCTGCCCGTGGAGAGGATGCACCTGGGAGGCAGCTCTAAGCC TTTACCTACCTCCTTCCCCAGGGGCTGCCTGCTGTGATGGTGGTGTGGAGTTTGGGCCCAATGTCACAG ACCCTCAAGATGTCACATCTAAGTGACCTGTGAAAGCAGCAGACAGACAGCCCGAACCCTATCTAG ACGCCTAAGAAAGCTGGATGGGTAATGCTAGTATGATTTCCACTTTACAACAGACGCTGAGGTAGAGA GCATTTGTTCTTATTCACAGCCAAGAAAAATACCAATTATTTCCAAATAAAGCAAAAAATTGGAACAGA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_020798.4
Summary:	This gene encodes a member of the peptidase C19 family of ubiquitin-specific proteases. These deubiquitinating enzymes (DUBs) catalyze the removal of ubiquitin proteins from other proteins. The encoded protein associates with polarized mitochondria and has been shown to inhibit NF-kappa B activation and delay PARK2-mediated degradation of mitochondria. Expression of this gene is upregulated by the let-7a microRNA and reduced expression has been observed in human tumor tissues. [provided by RefSeq, Jul 2017]
Locus ID:	57558

MW: 53.1