

Product datasheet for **SC210973**

URB2 (NM_014777) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	URB2
Synonyms:	KIAA0133; NET10; NPA2
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_014777
Insert Size:	933 bp

Insert Sequence: >SC210973 3'UTR clone of NM_014777
The sequence shown below is from the reference sequence of NM_014777. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CATGAAGGAGAGAAAAGATATACGGCCTAAGGCTATGGGACAGAAGTGCCGCCAGTGACACTGTCCAGA
GGCTTTGGCTGCATGGTCTGAAAGAGCTGGAGAATGAAAGACTTAAGATGTTCTAATTCGTAGTATTGG
TATACATAGAAAATCCTTTGGGGTTTATGTAGTATATTTTGATGTATTTTACATCGTGTCTTTCTTACT
ATTTTTTAATACATAGTTTTATGCAGTAAGTATTGCAATAGAATCCTGAAAATTGACCTGGGATGAGA
TTAATTCATAGAAAATTGCTGACTCTTGGGACCTTTCTGTGTTGGTTCTCGTCTTGGCTCAGTGGT
TTGGTGTCCCTCGTCTGCACTGGAACTACATAAACTTGGCTTTTTACTTTGGGTACATGGGCGTAT
AATTCAGCCCTGTTAAATATACTTGCTTTCAAATTTCTCAAGTAACATGGGAAGTATTCTTGAAATG
TCACATTTTCTGCCTTCCCTCTAAGTATGCTTTCTGAAGAAGTCAGGGAAGTTAGAGTCTGTGGCCTG
AGGTGTCTGCTCTGGGTGGCGATAGTGGGCACCTCAGGCAGGTCGGTGACGTTTAGCACAGGTGCCAGG
GCTCCTGCCTGCTCCTCTGTGTAGCTCTGTGAAGTTCATTTAGGAATTTTTTTTCTATGCAGTTT
AAGAAATAATCCTAATTGTTTTTTCTATTACCTAAGCAATATATTTTATTATAGCAACCTCAGAAAA
GAAAAATAAAGGATAATTTAAAAAATCATTATAGTCTCAGTTACCCAGATAACCTCGGTTGTCAAC
TTGGAGTATCTTGTGTAGTCCCTTACTATGTGTATGTATATAGATGTGCATATAAATATATATAGTA
GCTAAATTGGATCATAAATGCATTTTTTTAAAGTTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_014777.4</u>
Locus ID:	9816
MW:	35.8