

Product datasheet for **SC210977**

DDX26B (INTS6L) (NM_182540) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: DDX26B

Synonyms: DDX26B

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_182540

Insert Size: 904 bp

Insert Sequence: >SC210977 3'UTR clone of NM_182540
The sequence shown below is from the reference sequence of NM_182540. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AGTCACATCAACAGCAGATCATCATGTAGTGCAAAGACCAGTGAGAAAAAATGACAAGTTTTCTGTG
CTGTAGGATGGAACAGGATATTGTTGAAGCCTCCTGGAATGTTGAGTCAAGGGAATTGCTTCCAGAT
GCTAAGAAGCAGCAGTGGGGCTTTTGAATTTTATGATTATCTGGCAGTGAAAGCTGGGCTTTTGCCTTA
ATAATTTTTTAAAGTATGAATTGTTTTGTTTTGTTTTCTCAATTGAGGAAGCTGATGTTATTAATTC
CAGGCTAAATTCGGTAAACACCACTGCCCTACCACGGGTAATGAGAGGTCACTCACTTGAACCTTGCC
ATTCCAGGCATTCTCAGAGTGGCGAGGGGCCACCTGCAAGTGGAGCACAACTTGGTGCTCTTACTGTGT
CCTTCAGAAAGAATAGGTGTACAGAAAGGAAATGGCAATCTTATGTGTGCTGAACAAAGTTTTCAACAA
TTCCTAGTTGTGCCTTTTAAACCATGCAATATTCAAGGATAGTTGAATCAAAGAAGTAAGAAGCTGCTA
TTTGGGTAACCTTATTTCTGTGGGAAGGGGCAGGGAGAGTCACCAACAATCTACCTCAACTCTCTT
CTCTTTTGTCTAGAGACATTACAAAGTGCACTTGAGGCTGCCCCAACCTCTGACATTTGTTCTTGCA
GTGATGATAGAAAGTCTTCAGATGGACTTATACATTCTGTGCTTTGGAAGCACAGAAGAACAATAAT
GTGTATATTTCTTTAATGTTTATACAAAAGTTTATATGGAGCAGTATTGTTATGTTTGTATGAATTTG
CAAAAATTAAGTGTACAAAGAGATTTTGATTTGCATATATAAAATAAATCATTTTATTGATTTTCAC
AAGTTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_182540.7</u>
Locus ID:	203522
MW:	34.5