

Product datasheet for **SC210980**

NPRL3 (NM_001077350) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NPRL3 (NM_001077350) Human 3' UTR Clone
Symbol:	NPRL3
Synonyms:	C16orf35; CGTHBA; FFEVF3; HS-40; MARE; NPR3; RMD11
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001077350
Insert Size:	1349 bp



[View online »](#)

Insert Sequence:	>SC210980 3'UTR clone of NM_001077350 The sequence shown below is from the reference sequence of NM_001077350. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC ATTGCCGTCTTCCAGGCTCTGCTCCCCAGGCCAGGCGGAGGGCGGAAGGCTGCTGGGGTGCGCAGGT GGGCGCTCGCGTCTCCCCACCCAGGGCTCCCCCGTGCTGAGGCTGAGCCCTCTTGCCCTGAGGCCT GGCATTGGGTGGATGCGGGCTGGCCGTGGCCAGTGAAGCCTGCAGAGCCCCGCTGTCTTGCCCTGG TGGTTCGCTGTGGGGCTGCTGCCCTCTGTGTTCTACGCTTCCCTCCAGTCTTGCCGCACGCGTAGG CATCTCCACGCTGGAGGGGGCCTGCCAGGACTCTGTCTCTGGGTGAGGCCGCATCTTCAAGGCCCA TGTGGGCTGTGCGTTTCTCAGACCCTCCTTCTGTCCCTACACCTGCTCCTTGACCCCCAGTCTGTGG CCACCCTGAAGAATGTGCAGAAACACTTGTGTGGCCTGTCCCTGTCTCTGACAGCCTTCTATTTGTG AAGTGCCCTGTGGCCCCCTCCCCAGCACCTCTGTCTGCCATGCGCTTCTCCTCCAGGCTACCCCTGAG CCTTCTCCTGGCCAGTCTCACCACAGTCCACAGAAGCCGCAAACAGGCTCATAGCCAAGCTGTGACCT GGTCTGACCATCTGGGGCAGGAGGGCCTGGGCTGGCCCTGCTAGGCTGGAGAAGCCCTGTACCTGTG CACATCTTGCTGGTGGAGGCATGGCCACTGTGCAGGACCCACCTCGGGGGCTTTGCGTCCCACAC TGATGATTCTCCCAGCATCCACACGGGCCTGGCGCTACGTACTAGGCCCCAGCTCTCGTGTCTG GAGGAAGAGCTAGCTCCAGACATGGGTTGATCACCTAGAGGAGCTCTGGCTAAGGCACAGTTTCTAGA AATAAACATTTATTCGGTTTTTAAACAATTAATAATCTTGTCTGGGCAGGCCCTTTGCTCAGGCCTG TGTGCTCTGCTCAGCCACTCTGTGACCACTGACCCAGGGGCTGCCCCGACATAGAAGCGGAGGGG TTTCCGGGCCCCACTCCCCTGCATGGCCGACGCCACCCGGGCTGCTGCCACTACAGCCGGCTCACTGGG CTCCAGGGGACACGCTCCAGCCATACAGCTTCACTCTGTGCCAGGTCCCTCTGGTCAAAGCTCTTGTT GATGGCCAGGGCCTGGCAGCTTGGAGGGGCCACTGCAGAGCTCGCGGTCTTGAGGACACGGCTGGC GGTGCCTTTCCGGAGGGTGTGCGAAGCTGACGCATGGTCTCCAGACCTCCAGGGGCTCCAGTGTCTG CAGCAAGACGCAAGCCCCGTCCCTGTGGGACGGACAG AGCGGACCGACTTACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCC CAACCTGCCATCACGAGATTTGATTCCACCGCCGC
Restriction Sites:	Sgfl-RsrII
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001077350.3</u>
Summary:	The function of the encoded protein is not known. [provided by RefSeq, Aug 2011]
Locus ID:	8131
MW:	47