

Product datasheet for **SC210983**

C3orf33 (NM_173657) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: C3orf33

Synonyms: AC3-33

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_173657

Insert Size: 930 bp

Insert Sequence: >SC210983 3'UTR clone of NM_173657
The sequence shown below is from the reference sequence of NM_173657. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AGTCGCATAAACTTTTCGTAGAAAAGGGTGAACGGTCAAGAAATCTAGAGGTAACCTTCTCCAAAGAGT
GAAATATGTAATATATGGACGTTTTTCATTGAATCAAAATGGAAGTTTTGCCAAATGATCAAAGTGTT
AGTCCAATGGTATTTAAATATTTAAAGAGATGATTATTATAAGCATAATGTCAGAATTAATTTAAACA
AATTATAATTAATGTAATATGATTTTAGGAAGGATGAAACACTTTATGAGAAGTTACCAAGGCTGAATGA
TTTACATATATAGTGCTAATGTGTGAAGAGAAAGTTTAACTCTTCACCCACTATGTTGTTTTTAGTTAT
TTTGGGCTGCCATATTGAATTTTTGTTGTTTGTGTTGAGAAGGCGTTTCGCTCTTGTGCCCAGGCT
GGAGTGCAATGGCGTGATCTTGCTCACTGCAACCTCCTGCTTCTGGGTTCAAGCGATTCTCCTGCCT
CAGCCTCCTGAGTAGCTGGGATTACAGGCATGCACCACCATGCCCGGCTAATTTTGCAATTTTAGTAGA
GACGGGGTTTCTCCATGTTGGTCAGGCTGGTCTCAAACCTCTTGACTTCAGATGATCTACCCACCTTGGC
ATCCTAAAGTGCTGGGATTACAGGCGTGAGCTACCAAGCCCGGCTTTTTTTGTTTTGACATGGAGTTTC
GCTCTTGTCACCCAGGCAGGAGTGCAAGTGGCGTTATCTTGGTTCACTGCAACCTCCGCCTCCTGGATT
AAGTGATTCTCCTGCCTGAGCTCCCAAGTACCTGGGATTACAGGAGCTCCACCAAGCCTGGCTAAGT
TTTGATTGCCATATTGAATTTTTGTTTAGAAAAGATATGTGTTTATGTTCCAGTTTATTACCTTTT
ACTATTTATAATAAAGGAATCTTTATGCTTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_173657.3</u>
Summary:	Isoform 2: Secreted protein may play a role in transcription regulation via the MAPK3/MAPK1 pathway through an unidentified receptor on the plasma membrane.[UniProtKB/Swiss-Prot Function]
Locus ID:	285315
MW:	36.2