

Product datasheet for **SC210989**

RNF26 (NM_032015) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: RNF26 (NM_032015) Human 3' UTR Clone
Symbol: RNF26
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_032015
Insert Size: 915 bp
Insert Sequence: >SC210989 3'UTR clone of NM_032015
The sequence shown below is from the reference sequence of NM_032015. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ATCCTGCAGACCTCAATGTCTACCTCTAAGCCTCCTTCCCTGCCTGCCACCCCTCCATGCTCCACG
CAGGCACTCACGCTAGGACAGCATTAAACACCTCATCTCCGGGTCCTGGTCTGAATCCCCTCCTACCCCT
GTGGCCATCCTGCCATACATCCAGGACATTGAGTTGGAAGACTATGATCTGGGTGGGGCAGGATAACA
TGGCTTCTCTTTACCCAGTGGGTCCCTTCGATGCTGAGGGTGGTGAGTATGCTACTATGCAAGGGCCCT
GAGACTATTTGCTGTGGGCTCTCCTCCAGCCTGCCAGGGCCACCCAGATGCCTCTGGGGTTACCCCT
GTCTGCTCTGGTTTTTCTGTTGGAGATCTATAGGTCCTTTTCTGCCTCCTTCACATTTCTCCCCAG
CTTTTGGCGCCACAACACATCAGTGTCTTTGGGTGTTTTGGCAACTCAGGGGCCTTCGGATGATCTTA
AACCTTTGTGTTAGCCAGAGCCCTGTGCCCTGGTAGGCGTTGGGGTTAGTATCTCTCGGGTGCCCTC
AGAGCCACCTCTGCCTGTGATCGTCTGATGAGGCTCCCTCCCAACCTGATCCAAAGCCAGTCTCAGGA
GTTTACCCCTGGGATGGGGATGCATCTGCACCTGACTTTGGGGCCACGTGCCCTGTGGCACCCAGCT
CACTGGGAGTCTCAGGAGGGATAACCGGATTTCTGCTCTTCCCTGTCACTCCACATCACACAGAAA
AATGGCATTCTCTGTCTCTCCCTGGCATGGAGAGGGCAGACTGTGCACATTTCACTAGGGTCCAAA
TACAGAAAGGGCCAGGGCCAGGGGCTTGACGTTTCGTGAGGGGTCTCTGGCCAGTTTCCAATGAATA
AAGTTCTCTTGACAGCTC
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_032015.5
Summary:	The protein encoded by this intronless gene contains a C3HC5 type of RING finger, a motif known to be involved in protein-DNA and protein-protein interactions. The expression of this gene was found to be upregulated in cancer cell lines derived from different types of cancer. [provided by RefSeq, Jul 2008]
Locus ID:	79102
MW:	34