

Product datasheet for **SC211000**

DNAJB2 (NM_001039550) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DNAJB2 (NM_001039550) Human 3' UTR Clone
Symbol:	DNAJB2
Synonyms:	CMT2T; DSMA5; HSJ-1; HSJ1; HSPF3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001039550
Insert Size:	928 bp
Insert Sequence:	>SC211000 3'UTR clone of NM_001039550 The sequence shown below is from the reference sequence of NM_001039550. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GCTGGGAAGAAACCCGAGATGTGTTCTGAGCTGGATGCCGGTTCCAGAATCGTGCACAGTTCCAAC
AGGACAGCGCCTTCCCCATGCGCTGGGAGGGGACCTCCATTTCTCCCCTCACCCATGCTGAGTGTA
GAGCCGGGGCCTGGGTGGCGGTGGGGCCGGTGGGAGGTGGCAGTAGTCTTAGCCTGTGCACTCTCT
TCCTTGGGTGTTTGGTGCTGGCTCCTGGGACTACAAATCCCAGAGTGCAGTGTGCCCGGCCTATTTC
TGATAGATCCCGCTTGGGGAGGTGGTGTATGGTTACGGAGCTGTGCATCTTGGGACATGTAGTAGCCC
AGGTCCGCTTGTCACTCGCTGTGAGATGGGGAGATTTTGTCTTTTGATTATCCCTGTAGGGCTGGCAG
GGTTGTAGATGAAGGGGAATGATCTGAGCCTTGTTCCCTGACACGTCTTGCTAGCCCCAGGGTTAG
AGTGGGCAGGGCAGAGCCGCGCAGCACCTGGGAGCGGTACCTTTCCCTTGGGCAGCTGGGGTCCCAGG
AACAGCCAGGGCGAGTGGCATGTCTGCCTGAGCAGGTGTGGCCCCAGAAAGCTGAGGAGTGTGGGCT
GGCAGAGAGCTTCGAGGGCAAGGCCACCCGCGGGGCGTGTGTGGTGGGGCTTGGCATGTGATGGCA
GCTCCAGGCATGCCGCTGCTTGTATGGCTTTCTTTGGCCTCTGACCCTGCTGCCATTCTTTCCAACAT
CACAGATGAAGTGCCTCTCCTCCTCCCTGCCTGGGGAGCCAGTGCCAGGGAGGGGAGTGGTGGAGCC
AGTCGCTGTAACACTGAGCCTCAGAGACGAACCAACAGCTGGGCTGAGCTCAGATCCAGGGGGAAA
TGCTGGAAGTCAATAAACTGAGTTTTGAGA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001039550.2
Summary:	This gene is almost exclusively expressed in the brain, mainly in the neuronal layers. It encodes a protein that shows sequence similarity to bacterial DnaJ protein and the yeast homologs. In bacteria, this protein is implicated in protein folding and protein complex dissociation. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2011]
Locus ID:	3300
MW:	33.1