

Product datasheet for SC211036

BNIPL (NM_138278) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	BNIPL (NM_138278) Human 3' UTR Clone
Symbol:	BNIPL
Synonyms:	BNIP-S; BNIP-Salpha; BNIP-Sbeta; BNIPL-1; BNIPL-2; BNIPL1; BNIPL2; BNIPS; PP753
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_138278
Insert Size:	1063 bp
Insert Sequence:	>SC211036 3'UTR clone of NM_138278 The sequence shown below is from the reference sequence of NM_138278. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CGGGATCTCCATGGCTCAGGAGGGACAACACAGGACTGGATAAAGGCCTTAGAACAGTTAGTGAT
CTGCCTACACCTGAATCCCTGAAACATCTGAAGTGTGTTTGTAAATCATCTTATCCCCAACCTCAGTACC
ACCGGATCTTCACTTCTCAGTGGGATTTTGTCTTTGCATGACCCTACTTTCAGCAGGGCTTTGAGGCT
GGGAAACTGATCTGCTGGGTGACTTTTCATTCCAGTTGCTGGGACGGAAGGCTGAATGTAGGGTCATTTT
GTATGGGATATGCAGAGCTCTGGGTTTTTTGTTTTTTGTTTTTTGAGACGGAGTCTCGCTCTGTCGCC
AGGCTGGAGTGAGTGGTGCGATCTCGGCTCACTGCAACCTCCGCTCCCGGTTCAAGCCATTCTCTG
CCTCAGCTCAGCGAGTAGCTGAGACTACCGGCGCACGCCACCACGCCAGCTAATTGTTTTGTATTTT
TAGTAGAGACGGGGTTTACCATGTTGGCCAGGATGGTCTCGATCTTTGACCTCGTGATCTGCCCGCC
TCGGCCTCCCAAAGTGCCGAGATTACAGGTGTGAGCCACCGCGCTCGGCCAGAGCTCTGGGTTTTAATC
CTACTTTAGCTGTAGAACCTTGGGCAATAACTTCATCTTTTTGGGTCTAATTTTCTCCTGTCTAAC
TGGAGGGACTGTGATCCTTCCACCTTGGATGGAAGAGGCTTACCTGACAGCCAGCCTGCCTGCTGGGAG
CCAGGAGAGTCAGCTCATTAAATCTTGAAGAACCCTGTATGCCCTTTTCTTCTAAGTCATGTCTGC
TGCTGTGAGCCTGGGAAGGAGTGCTTTCAAACCTGTATTTTGGCCTCTTTAGTAAATTTAGAAGCTG
TAGAGGCTAATAAACTGCGATGAGACATTTAAGCTCTGCTCCAAAGAAGTGAACGACTCTTTCACCACC
CCTTGATCCCAAACCTTCGGTTCTGGCAGAGTTCTTGCCGAGAGGTTCTTAGGAGCACCCCGCGCG
GCCCGCGGAGCGCGCTGCGCTCGAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCGCCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_138278.4
Summary:	The protein encoded by this gene interacts with several other proteins, such as BCL2, ARHGAP1, MIF and GFER. It may function as a bridge molecule between BCL2 and ARHGAP1/CDC42 in promoting cell death. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2011]
Locus ID:	149428
MW:	39.3