

Product datasheet for SC211052

XKR8 (NM_018053) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	XKR8 (NM_018053) Human 3' UTR Clone
Symbol:	XKR8
Synonyms:	hXkr8; XRG8
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018053
Insert Size:	923 bp
Insert Sequence:	>SC211052 3'UTR clone of NM_018053 The sequence shown below is from the reference sequence of NM_018053. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GATGAGGCTGCTTCGCCAGTGAAGGGAAGGTGAACGGCGTCCTTTGAAGCAGGATCAGACCCAGCCAG
CAGAGATGGAGAGTGACTCTGTTGGCAGAAGGCAGGCGAGGATAAGCTAACGATGCTGCTGTGCCCTCT
ATGCACTCAGCAAGAGCGGGACGCTGTGCTGGCCGGGCACCAGGGATGGTGTGAGTCGGGCAGAGG
CCTCCTTTCAAGGAGTTCACAGTGAACAAGATGAGAAGGGCTGGGCCCTGGAGGGTCAAGAGCCCCAAT
TATGTACAAGACACTTTGGGAGGAAAGAAGACTACCTTTTCCCCCTGCCATTGGTATAGCTGGTGCCCC
AAAACCTCCACCTCCCTCCCTGGCTACCTCTAAAATGACTGGTATAGGTGCTGCCCCACCCCTTAGCTC
CCCTATCTGGGCTAGGAGGCCACAGGGGCTGTCTCTAGAATTCTTCTTCCCTCCCCCACACCATTC
ATTCAATTCATGAAACAAATCTTTGCCAAGAGCAGTTTATGTGCCAGGAACATCATTCTGTCTTGCAA
CCTGGAACAAGACCAGCTACCAGCCTAGCTTCATCCGCTACTTGACCAACACAGTCCCGGGTTAGATCC
CAAATGCTAGAAGCCAGGGATGCCCACTCTGGGTGGCCCCAGTCAGAACCTCTGGGATCTCAGTGAAG
CTGGCCTGGCCTCTGCTCCTGCTCTCAAGGGGCTGCTTTTCAACCAAGAGCCTTTGTGAGCCTGGTCTGA
GCCTTGACAGCCACTGAGTATTTTTTGCCTTAGCCAGTGTACCTCTACCTCAGTCTATGTGAGAG
GAAGAGAATGTGTGTGCTGTGGGTCTCTACAAGTGACAGATGTGTTGTTTCAACAGTATTATTAGGT
TATGAATAAAGCCTCATGAAATCCTC
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_018053.4
Summary:	Promotes phosphatidylserine exposure on apoptotic cell surface, possibly by mediating phospholipid scrambling. Phosphatidylserine is a specific marker only present at the surface of apoptotic cells and acts as a specific signal for engulfment. Has no effect on calcium-induced exposure of phosphatidylserine. Activated upon caspase cleavage, suggesting that it does not act prior the onset of apoptosis.[UniProtKB/Swiss-Prot Function]
Locus ID:	55113
MW:	34