

Product datasheet for **SC211069**

CENPL (NM_001127181) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CENPL (NM_001127181) Human 3' UTR Clone
Symbol:	CENPL
Synonyms:	C1orf155; CENP-L; dj383j4.3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001127181
Insert Size:	928 bp
Insert Sequence:	>SC211069 3'UTR clone of NM_001127181 The sequence shown below is from the reference sequence of NM_001127181. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ACAGAAGTGGCAATTTTCAAATTGAGTAAGCCTTATGTGGACTATAAGTTATAGATTATATACTCTT
ATTGATAACTTGCCTAATTGCTATGCTGAAAGAGACTGCAGGAGAAATAGGCATCTATCTCTGCATCTG
TTTTCCCACCATGCCTTTGGAGTTGCCAAGATGGAAGCCAAGAAGGATCTAGAAGAACAAGAATATG
GTAGTAGATGAGCCACAGCCAGGTGCCCATGTACTAATCATGATAACCTGACATGCCATTCTCAAAATG
CTGAGTTGTTAATTTCTTGTCATCTTTAAATATATATATATAGGCTGGGCTTGGTGGCTCACACCTGTA
ATTCCAGCACTTTGGGAGGCTGAGGTGGGTGGATCATTTGAGGCCAGGAATTCAAGACCAGCCTGGCCA
ACATGGTGAAACCCCTTCTACTGAAAATACAATAATTAGCTGGGCGTGGTGGCACATGCCTATGATC
CCAGCTACTGGGAGGCTGAGGCAGGAGAATCGTTTTAACCAGAGAAGACAGGCTGCAGTGAGCCAAGAC
TGCACCACTGCACTCCAGCCTGGGCAACAAAGTGAGACTCTGCCTCAAAAAATAAAAAAGAAATAAAAT
AAATAAATACAGCGTTTATTATGTTTAACTCTTAGGATTCAAATGAGATGTATCTAATGGATATTTCT
GATTTGGGTCTTTAAACCTTTAATCTCTTTAAGCATTTTAGACTTGCATTAGAAAGAATTCTGTCT
TCTAAAGACATTTTCTAGCATTTCACTGCGGAGAACTGGCAATCAAAATCCCCTCAAAAGGGAATAA
AGATTAATATACACATTTTATAGTCATTAATATGCTTCTGTTTTAAATACAGTTCAACTCCTGTTGCA
ATGTTGAATTAATGTTTATTTTAAATGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001127181.3</u>
Summary:	CENPL is a subunit of a CENPH (MIM 605607)-CENPI (MIM 300065)-associated centromeric complex that targets CENPA (MIM 117139) to centromeres and is required for proper kinetochore function and mitotic progression (Okada et al., 2006) [PubMed 16622420]. [supplied by OMIM, Mar 2008]
Locus ID:	91687
MW:	35.4