

Product datasheet for SC211072

C6orf192 (SLC18B1) (NM_052831) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C6orf192 (SLC18B1) (NM_052831) Human 3' UTR Clone
Symbol:	C6orf192
Synonyms:	C6orf192; dj55C23.6
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_052831
Insert Size:	930 bp
Insert Sequence:	>SC211072 3'UTR clone of NM_052831 The sequence shown below is from the reference sequence of NM_052831. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CGAACTACTCTCTTGCCTAATGAAACCTAGTCCGATGGATCCTGGATTGATACAAGGTTGAGAAATGAA
TGCTCCTGGCCTTAAACATCACCGTAGGAAGGGTTTTTAAATTTTACGCGCAAACTCCGTGGACCCC
GTGCCAGTGTCTTGAAGTGTCAACGTGTTTTTGGATGATCCTGTATTGGGCTGTACTTACTGTGATAC
TGAAAAGCTGCTGCTGAAGCAGCTATATTTGAAATATTAAGTATGAAAGGAGTAATTAACAAGC
AAAACAAAACAAGACTTAGTTTTTAAATGACCAAACTTGCTTAAAGATGTTGTTATTAACCTCGAGTT
AGTTCTTATTTCTCTGTTTATTTTTATTCTAAGTACACTGATTCTGTGAATGTACCTTTTTTATTAA
CAGGGAAGAAATGAATTAATTTGATATGCTCTAAATACATAAAGGTGCTTCAAATATGTAGAAACAT
TACTATGAAATCAGTTTTTAAAGATATACTTTCTTTTGTCTGAGGTTTTTCGGTCTTGTTCAAAAG
GAAGAATTCTTGCTGCCATACAGAACTCTTAGCACTCCCTGACCTTAAGCTTTTCTAAAAATTCTG
TTTGTGTGAAAAGTACAAGAATAACAATACTTACAACCTCCATTTTGTAACTACGTTCACTTATGAT
CTGGATTATATAACATTACTTGGTATAACGTTTTTCATTTCTTTAATGTCTCTGTTTTTGGCTCTAC
CATCTGTTTTGTTTTGTTTTATCTATATCTTGGTAGATGATTTTCATCCCTAGAGCAGGTCAGCCTC
CTTCCCCTAATGCGAATGCTTGTGTTTGTAGGGAAGGGCTTCTCCAACCTTCGTGTGAAATTGTGATGT
TGAAGTGAATAAATGTCTATTGTGTAACCTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_052831.3
Summary:	This gene encodes a protein, which has high sequence similarity to rat, xenopus and zebrafish proteins. The protein function is unknown. [provided by RefSeq, Jul 2008]
Locus ID:	116843
MW:	36.5