

Product datasheet for SC211096

MRGPRF (NM_145015) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MRGPRF (NM_145015) Human 3' UTR Clone
Symbol:	MRGPRF
Synonyms:	GPR140; GPR168; MRGF; RTA
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_145015
Insert Size:	910 bp
Insert Sequence:	>SC211096 3'UTR clone of NM_145015 The sequence shown below is from the reference sequence of NM_145015. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATGCAGTGTCCCCGGGAACGCCTCTGAGACTCCAGCGCCTGGAGGAGGCAGGGGAGGAAGCGGCC
TCCAAGACCCCTTCGCCTTGGGACAGGAATGGGCACCTGCTTCTGAGTCCATACAGGAGAAGAAAGATCT
GTTTCCTCTCCTCGGCCTCCTTCTCCCTGGGCTGGGACTCCAGGGTGGCTGGGAGACTGGGCAGCC
ACCAGCAAACAGACCCTGTGGCCCTGCCCGCTCCCCACCCATTCTGCTCCCCTAGAGACCTCTTGT
ACAGAAGTTGCCCCAGGTGGTGGGCCCCCTCCTTGCCCTAGGCTGGTTGGTAAAGAGAGAGGTCAA
CACCCAGCCTAGCCACCTCTGCCTCTTGGGTACGCCCTCCTTGACTGTGTCCAGCCAGCACCAGGCCA
GCAGCCTCATCCCTGCCATTGAGGCTGTTCCAGAGATTCGATCCTCTTAAGGCATTATCAGTGAGCAA
ATGTGAAGGAAATGGTGTCTGGAAGAAAGTTCTGGTTCACATGCCTTGTAGCTAAGTCTTTTGCAAAC
AACCTCCCTTCCCCCGTCGAGTCATTTGGTGACTTTGATGGGGGATTCTGTTATGTCAAGGCTCT
GGAGACAGGAAGGGCCTTGGCCGCCTTGGGTAGTTGACCTGCCTTTTCTGACTCCGGGACGAGCCAGT
CCTAGGCTGCCTCCGGGAGCACTTGAGGTATCCCGCAGGCCATGAGGACCACTGGGCAGCTCCTGGAC
AGCCTCTTGGCTCCAGCCCCACCCGAAAGTGGACACTGGCTCCGCCCTGGCCACCTGGGACTGGCAC
TGTGGTGACAGTGCCCAATGTGGCCAACGGAAGTTTTATAAAGACAAATGTATATCAATAAACAT
TTTATAACTTGCA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_145015.5
Summary:	Orphan receptor. May bind to a neuropeptide and may regulate nociceptor function and/or development, including the sensation or modulation of pain (By similarity).[UniProtKB/Swiss-Prot Function]
Locus ID:	116535
MW:	33