

Product datasheet for SC211099

ATP9B (NM_198531) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ATP9B (NM_198531) Human 3' UTR Clone
Symbol:	ATP9B
Synonyms:	ATPASEP; ATP1IB; hMMR1; HUSSY-20; NEO1L
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_198531
Insert Size:	930 bp
Insert Sequence:	>SC211099 3'UTR clone of NM_198531 The sequence shown below is from the reference sequence of NM_198531. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CCTCCCAGCTACTGCAAGCTGGCCTCCAGGGGCTGTGCACCCCCAGCGGGCTGGCCCCAGCACCTTC
TGCCCTTCCCAGCACCTTGTGCCCTTGCCAGTGAACGCAGGGTTTGCCATTGCTACCAAGCAAGCACCA
CAAGAAAGGGAGGGTACGCCAGGCGAGCCAGGGCACAGATGCTGAGACAGCCTCTCCTTCTCAGTGCA
GGGACGTCACCCCTGCCAGGCAAGCCCAGGGCACAGATGCCAGGATGGCTTCTCCCTCTCAGTGCGAGG
CTTCAACCCTGCCAGGCAAGCCCAGGGCATAGATGCTGAGACAGCCTCTCCCTCTCAGTGCAAGGACGT
CACCCCTGCCAGGCAAGCCCAGGGCACAGAGGCCGGGACGGCCTCTCCCTCTCAGTGTGAGGCTTACC
CATGCTAGGCAAGCCCAGGGCACAGATGCCGGGATGGCCCTCCCTCTCAGTGCGGGAACGTCACCCCT
GCCAGGCAAGCCCAGGGCACAGATGCTGCGATGGCCTCTTCTCTTAAGTGTGGGGCCTCACCCCTGCT
TTTCTTCTTTTTTTGTATTGTCAAAATTGATTTCCATATTGAAGCAGCTTGAGTTTCTACTGAAAT
GAGCCCGAATTATTTCACTATTACTGTAAAGGGTTCATCTTACTCTGGCATTCTGAGAATTAGACTGAA
AGTTTAATTTCTGCAGTTCCTCATATTCAGATTCTTTCTTTGATGTTATAACACAAAGTCATTCCTAC
TCAAATGTAATAAAATTGAGGCTCCACGGAGGCCCGTGTCCACATCAGCTGCAGTCCCCTCAGACAG
ACGGCCTGCCCCACAGCCTGCCCTTTGCTGAAGATACCCTCGCTTTATACTTGCATTTGCTTTTCAGT
AAAAACAGTCAAGTTACTGAGTTCTCACATAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_198531.5
Locus ID:	374868
MW:	33.4