

Product datasheet for **SC211113**

UMODL1 (NM_001004416) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	UMODL1 (NM_001004416) Human 3' UTR Clone
Symbol:	UMODL1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001004416
Insert Size:	952 bp
Insert Sequence:	>SC211113 3'UTR clone of NM_001004416 The sequence shown below is from the reference sequence of NM_001004416. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AACTTCAGCTACCAGGTGTTCTACGAATAGGAGGCGCAGGCTGACAGGAAGGTCGCCGTGAGTCAAGCT
GCCTCCAGAACCTCAGAGCTTCCCTGGTGGGCTCCCCGGGATCCCCAGTGTCTCTGACACCTCCACC
CATCCCTCGGTTCTTAACTCTTCAAGCCTTAACGGAGGTCTGCTCTGACGGGTGGGCTCTGCCAGAGCC
CGGGTGAGCCCAGAAAGGAAGACAGCAGCCACCGTCTGTCCCGAAGAGGCAGGCCGTCTGTAGGTCCT
AGAGGAGCCACAGCCCAGGGGCAGATGAAGGGGCTGCGGAAGACGGGGGCAGTCCTGGGGGTGCTGCGG
CTACACCACCACCCGCGCGGCCCGCAGCCTAGACCTCCAGGCCTGTGACCCTCCACACCAGCCCTC
AGAACCCTCCTGGGCTTGCCCTCCCTTGCGTCCGTACCCTTTGGCAAATATAGAATATTTACATTC
TCAGAGAGACCCGACCGCGCTCTTGATGCTCTTTCGAAAATAGGTCAGTCTAGAAATACACTGCTAAT
GTTATTTTATGTTGATGTTTATGCTGTTTGACTTTTCTCTGTGTACCAAGGTATTGCTTTTATTACA
CGACAGCGACTCAAAAGGCACTCGATTAATGTGACAACCTTTTCAATAAGCAGAAATAACGTAGGTACA
CATCACTCTTTACATTTTCTAAGCATTTTACAGCCGTTTCTTCATATAATCCGACCACAGTGGGAGG
TGTGATTTACCCATTACACAATGAGAAACCAGAGGAGCCGATGAGTTACTTAATTGAGGTCACAGAAATG
AATTAGCAAGAAAATGGTTCTAAATCTAAGTATTTTAGTCTAGAATTTTCTCCATTACATCATCCTAA
GAGATAATGCTCTGTACTTCATTTGAAATAAACTGGAATTGTATTAGATAGCTCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001004416.3
Locus ID:	89766
MW:	35.1