

Product datasheet for SC211120

FAM89A (NM_198552) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FAM89A (NM_198552) Human 3' UTR Clone
Symbol:	FAM89A
Synonyms:	C1orf153
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_198552
Insert Size:	935 bp
Insert Sequence:	>SC211120 3'UTR clone of NM_198552 The sequence shown below is from the reference sequence of NM_198552. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGCAGCGACTGGATTCTGGAGTCCATCAGAGGGTCTTGGGAGGGATGTGACTGTTGGGAAGCCCTTCC
TACTGGACACGCTGTCATCATTTGCTGCTTCTCTTGCAAGAAAGCACCTCCGTTGTGGACGGTCCTCGG
GCACAGGGGATGAGCGCTACCAAGTTTCATTTGTAGGCAGGGAGTTCTCCGCGGATGCATGGTGGCAGTC
TGCTTTGATGGCAGCAGTTTCTGCTTAGGTGACCTAGAGGTCCTCAGCAGTATCCTCCACACCTATTTA
TTGAGGTGCACCTGCTGGGGATTCATAATGAGAAATATAACAAGAGGATCTCGGTGAAAGGCCTTAGTGG
GTGTTTTGTGTGAGGTGGCTTGTAGCTAGCTACTTCCTTACAGATGGTAGAGTATTCCAATCCTCTTTG
TGTTAGGGTTCTTGCTTCCAGTTTGGGATGTATTAACCAACCATTTCACTGCTTCCCTTCTCAATAT
GCTCTGCAGCTTTTCTTGCTGTTTAAACCTCTCGCCTCAGCTTTATTTATTTGTAAGCTGCATTACTAA
CTGCCCAGTGATTCGGTGAAAGCTTTTACTGAAAAAGTTAACATTTCTAGTCATCCCAATCAACTGGC
TTTTTCAACCAAAATTTATATCATTCTTTGTCTATCAGATACGAGAGGAAGGAAGATAATACGAAGA
CATGTTGAATAGTAAAAAAAAAAAAAAGAACACAAAACTGGGGCAAGCCAATGTGATGTATCACTCAC
TGTAAGATGGCAAATGTTTTTCAATTTTAAAGATCCGAATGTAACTAGTGTGCTAGAAAGCAAACCACC
CGCCACTCAAACCAGTAATTACCTTAAGCCTTAATATATTATTTAAAAATACCTTATGAGAACATTACAC
TTGTAGGTAAAAATGAGGATAAAATGCTAAACTATC
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_198552.3
Locus ID:	375061
MW:	35.7