

Product datasheet for SC211131

TCCR (IL27RA) (NM_004843) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TCCR (IL27RA) (NM_004843) Human 3' UTR Clone
Symbol:	TCCR
Synonyms:	CRL1; IL-27RA; IL27R; TCCR; WSX1; zcytor1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004843
Insert Size:	958 bp
Insert Sequence:	>SC211131 3'UTR clone of NM_004843 The sequence shown below is from the reference sequence of NM_004843. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GGGCCCCCAGGCCACAGGTTCTGGCTGACCACACGTCTGGCTGGGGGCTGCCAGCCAGGCTAGAGG
GATGCTCATGCAGTTGCACCCAGTCTGGATTAGCCCTCTTGATGGATGAAGACACTGAGGACTCAG
AGAGGCTGAGTCACTTACCTGAGGACACCCAGCCAGGAGCTGGGATTGAAGGACCCCTATAGAGAA
GGGCTTGGCCCCATGGGAAGACACGGATGGAAGGTGGAGCAAAGGAAAATACATGAAATTGAGAGTG
GCAGCTGCCTGCCAAAATCTGTTCCGCTGTAACAGAACTGAATTTGGACCCAGCACAGTGGCTCACGC
CTGTAATCCCAGCACTTTGGCAGGCCAAGGTGGAAGGATCACTTAGAGCTAGGAGTTTGAGACCAGCCT
GGGCAATATAGCAAGACCCCTCACTACAAAAATAAAACATCAAAAACAAAAACAATTAGCTGGGCATGA
TGGCACACACCTGTAGTCCGAGCCACTTGGGAGGCTGAGGTGGGAGGATCGGTTGAGCCAGGAGTTG
AAGCTGCAGGGACCTCTGATTGCACCACTGCACTCCAGGCTGGGTAACAGAATGAGACCTTATCTCAA
AATAAACAACTAATAAAAAAGCAAAAAAAAAAAAAAAAAAGAAAAGAAAAACACTGCATTTGGGCACCAT
CTCAGCTCCCTTGATCCAGGTGCAGCATGGACTGAGTTCTTGACAACAGAATGTGGTCAGAAGTGACA
TATGCCAACACGGGGTCTGGGTGGGGCTCCCCACATCCTTTCTTGCCATGAGCTGGAACATAACA
CATGCCTATGATCCAGCTTTGGTCATACCCAAGGGGAAGGTGGAGCAAGAAATGAAAAGGAACCTGAAT
CCCTGAATGACTGCATGGATAGAACCACTAAGAAAAATAAACTTTTATATTTTATAGCTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_004843.4
Summary:	In mice, CD4+ helper T-cells differentiate into type 1 (Th1) cells, which are critical for cell-mediated immunity, predominantly under the influence of IL12. Also, IL4 influences their differentiation into type 2 (Th2) cells, which are critical for most antibody responses. Mice deficient in these cytokines, their receptors, or associated transcription factors have impaired, but are not absent of, Th1 or Th2 immune responses. This gene encodes a protein which is similar to the mouse T-cell cytokine receptor Tccr at the amino acid level, and is predicted to be a glycosylated transmembrane protein. [provided by RefSeq, Jul 2008]
Locus ID:	9466
MW:	36.2