

Product datasheet for **SC211148**

TBCCD1 (NM_018138) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TBCCD1 (NM_018138) Human 3' UTR Clone
Symbol:	TBCCD1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018138
Insert Size:	940 bp
Insert Sequence:	>SC211148 3'UTR clone of NM_018138 The sequence shown below is from the reference sequence of NM_018138. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GCAGCAGGCTCCAAACAAGCAGCTGGATAAGGATCTTACATGCAGACACTGGGCCCTGGAAATACAAGG
ATAAATGGAGGATGGTGCTCATCTTGAACTTGAACTAAATTGGGAAGGTGCTATTGCTGCTTAAGAG
ACATTGGGACTTTTCTGTCTTTTTTTTTTTTTTTTTTTTGGAGCAGAGTCTTGTCTGTGCGCCAGGCTG
GAGTGCAGTGGCGGATCTCGGCTCACTGCAAGCTCTGCCTCCCGGGTTCACGCCATTCTCCTGCCTCA
GCCTCCTGAGTAGCTGGGACTACAGGCGCCGCCACCACGCCCGGCCAATTTTTTTGTATTTTAGTA
GAGACGGGGTTTCACCGTGTTAGTCAGGATGGTCTCGATTTCTGACCTCGTGATCCGCCCATCTCGGC
CTCCCAAAGTGCTGGGATTACAGGCGTGAGCCACCGTGCCAGCCTTCTGTTCTATTTAAGAACTTTGA
TTTCCATTCCATTTTCATTTATTACTGTTTAGGCCGAATTTTCGCTAACTTCTGACTTTCATATAAAC
ATATTGCTTATAGAACATCAGAAGACTCGTCTTTGTGATGCCAAATGAAAAAGTTGGTGAAGGTAGTA
CCATTGTATGTATGTTAATATTAAGCAACTGAGTGATAGATTTGCAACTATAGTGTATCTTCATTG
TGTGATGAAAGAGTTTTATTCTTTCTAGAGTTAAATAAGAGTTTTATAGTTCTTTGTAACATTTAAAA
GGAAATTGCTACTCCTTTTAAAAAATAATTGAAAGGAAGTATGATTTGTTGGTATATAGACATAGAAA
TAGTTTTATACTATATTTAGTCCTTAAAGTATGGAACCTAATGAAATTGTTAAAAAAATCTTTATTT
TTCACTTCCTGCTTCAAATTTTCAATAAACAACTCTAAAGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCGCCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_018138.5</u>
Summary:	Plays a role in the regulation of centrosome and Golgi apparatus positioning, with consequences on cell shape and cell migration.[UniProtKB/Swiss-Prot Function]
Locus ID:	55171
MW:	35.5