

Product datasheet for **SC211151**

ANKRD20A2 (NM_001012421) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ANKRD20A2 (NM_001012421) Human 3' UTR Clone
Symbol:	ANKRD20A2
Synonyms:	MGC176486
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001012421
Insert Size:	957 bp
Insert Sequence:	>SC211151 3'UTR clone of NM_001012421 The sequence shown below is from the reference sequence of NM_001012421. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TATTTTCAAACCTCTAGAAAGAAAATTAAACATTTGGTTCTGGATACATGTTGAACCTAGTTGAATAT
AAAAATCTAGATTAAGTGTGTTTACCATACTGTATAATTCCATTTACATGAAGCATCCAGAAAAGAT
AAATGTATAGGGACAAAAGTAGATTCTGTTTGAAGGGCTGGGGCTGGAAGCTGGTAGTGACTGCT
AATGGGCATGAGGAATCTTACAGTGATGGAATGCTCTAAAGTTGGATTGTAGAGATGGCTGCACAACT
CAGTAAATGTACTAAAAATCTTTAACTTAAACAGATACATTCTATAGTATGTAATTTATTTCAAC
AAAGCTGTTTTAATAAAAAAGGAAAAATGTGTTTACTATATCGGCTTAGAAACATGCCTCATTTCTAG
GAAATAAAAGATAGAGGTGAGAGATGATTTACTTTGAGAAAAGACATTGTGTACCTATGAAATTTTAT
TAGGCACAGAGTCATATTTTAAGGTAGATAGTTCTGTATTGCTGAAATAGTAATTTAATGTCTTTATG
TTGCCACATGTTAAGACCATAATGTAGTTATAAATGGAATGTTTACACCTGAAGTGAGTATTTCAAA
TTAAAATTTAATTAAGTGATTTTCTTCGACACTTAATTTAGATTCCCGAGATGAATTGAAGTGATTG
CTGTGCTTGTAAATACCTTGCTTTAACTAGCTTTTATGTATTTTAGTTGGTATAGCTTTGTATTATT
CATATTAACAAATCTGAAAATATGTCAAATTACGTGTTTTATGACCATGTAATGTTTAAAGGCACCT
ACTTGTTATAAAATCATATTTTAGGATAAATGTGGTAAACTTAGCAAACTATATTTGGTTTAGTTTT
CCCACTGGTATTTATAGTTTACTTTGAATATTTATATTAATAATTAGCTCATAATTTTTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001012421.1</u>
Locus ID:	441430
MW:	37.6