

Product datasheet for SC211170

CWF19L1 (NM_018294) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: CWF19L1

Synonyms: C19L1; hDrnl; SCAR17

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_018294

Insert Size: 959 bp

Insert Sequence: >SC211170 3'UTR clone of NM_018294
The sequence shown below is from the reference sequence of NM_018294. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GAGCCCTATGACTTTACTCTGGATGACTAAACAAAGGGAAGAACTTTTATGAACTCCACAGGAAGTA
GTAAAGCTTTTTTTTTTTTAAATTAAGAATTTTTTTGAGACAAAGTCTCGCTCTGTCACCCAAGC
AGGATTGCAGTGGCATAACTGTGGCTCACTGTAGCCTCAACCTCCTGGGCTCTAGAGTTCTCCACCT
CAGCCTCATGAGTAGCTGGGACCACAGGCGCATGCTACCATGCCTGGCAAACTTTTTGATTTTTTATA
GAGACAGGAGGGTCTCCCTGTGTTGCCAGGCTGGTCTGTAATGCCTAGGCTCAAGGGATCCTCTGCCT
TGGCTTCTTAACCTGCTGGGATTACAAGCATGAGACACCATTCCTGGCCTAGAAGCCTATTTTAAAGA
AACTACAATCTCCCATGGGACTGTTCCCTGCCTCTTTGTGAGTCCCATGGAAGTTGCCTACAGCA
AGAGGCCTAAGATTGAATCTTTTTGGGAAAAAGTCATTCTAGGATGAAAATCCTATGTTAAGGCCGGGC
GCAGTGGCTCACGCCTGTAATCCAGTACTTTGGGAAGCCGAGGCAGGTGGATCACCTGAGGTGAGGAG
TTTGAGACCAGCCTGGCCAACATGGTGAAACCCCGTCTTTACTAAAGCTACAAAAATAGCTGGGCGTG
GTGCCAGGCACTTGTAATCCAGCTACTCAGGAGGCTGAGGCAGGAGAATTGCTTGAGCCTGGGAGGTG
GAGGTTGCAGTGAGCCAAGATCGCTCCATTGCACTCCAGCCTGGGTGACAGTGAAACTCCATCTCAAAA
ATAAAGAATAAAGTATGTCTGTCTCATCCAGCTCCTATGTCTGTTATCCAGCTCCAAGTACAGCTTGTG
TATATCAACATTTCAAAAACCTTTAACTACCAATGAAATCCAGATGTTTTCCATGGTT
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_018294.6</u>
Summary:	This gene encodes a member of the CWF19 protein family. Mutations in this gene have been associated with autosomal recessive spinocerebellar ataxia-17 and mild cognitive disability. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2014]
Locus ID:	55280
MW:	35.7