

Product datasheet for SC211203

PLAC8 (NM_016619) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PLAC8 (NM_016619) Human 3' UTR Clone
Symbol:	PLAC8
Synonyms:	C15; DGIC; onzin; PNAS-144
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_016619
Insert Size:	963 bp
Insert Sequence:	>SC211203 3'UTR clone of NM_016619 The sequence shown below is from the reference sequence of NM_016619. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AACAGAAGGAGAGCCATGCGTACTTTCTAACTGATGGTGAAAAGCTCTTACCGAAGCAACAAAATT
CAGCAGACACCTTTCAGCTTGAGTTCTTACCATCTTTGCAACTGAAATATGATGGATATGCTTAAG
TACAACTGATGGCATGAAAAAATCAAATTTTGGATTTATTATAAATGAATGTTGTCCCTGAACCTAGC
TAAATGGTGCAACTTAGTTTCTCCTTGCTTTCATATTATCGAATTCGAATTCCTGGCTTATAAACTTT
TTAAATTACATTTGAAATATAAACCAAAATGAAATATTTTACTGATAAGATTCTTCATGCTTCTTTGCTC
TCCTTAAATGTCTTTTCACTAGTTAGTTCCAAGGGTCAGTCTCATAATTTTGTCTTATACTTTGAT
TTCTTTTTTCTTTTTTTTTTTTTTAATAGAGTTGGGGGATCTTCTTTCTTTTTTTTTTTTTTTTTT
GAGATAGAGTCTCGCTCTGTTGCCAGGCTGGAGTGCAGTGGCGTGATCTCAGCTCACTGCAACCTCCA
CCTCTGGATACAAGCGATTCTGCTGCCTCAGTCTCCGAGTAGCTGGGATTACAGGTTCCCGCCACCA
TGCCAGCTAATTTTGTATTTTAGTAGAGACGGGTTTCACCATCTTGGTCAGGCTGGTCTCGAACT
CCTGACCTCGTGATCCACCTGCCTCGGCATCACAAAATGCTGGGATTACAGGCATGAGCCACCGCTCC
AGCCTAGAGTTGGGGTCTTTCTATGTTGCCAGGCTGGTCTTGAACCTCCGGGCTCTAGCAAGCCTCT
CATCTCAGCCTCCCAAGGTGCTGGAATTACAGGTGTGAGCCACTGCACCATACTTTGGTTTCTTAGCAG
GAGATTTTGGTTACAGTAACCCAGCTCTCTTTCATCTAAGAACTTTAGATGAGAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_016619.3</u>
Locus ID:	51316
MW:	37.9