

Product datasheet for **SC211218**

C14orf106 (MIS18BP1) (NM_018353) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C14orf106 (MIS18BP1) (NM_018353) Human 3' UTR Clone
Symbol:	C14orf106
Synonyms:	C14orf106; HSA242977; KNL2; M18BP1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018353
Insert Size:	949 bp
Insert Sequence:	>SC211218 3'UTR clone of NM_018353 The sequence shown below is from the reference sequence of NM_018353. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TATTATTTTTCGAACTCTGATTCTGCAATGAGAGAGAAAATATGATTCTGGGATTTTTACCATAA
AGCAGACAGTGTGTATTTTCACTGGAGTACATGTATTTCTTTGTAAGTAGCTTCCTATGAAAAT
GTGGACTTTTTGAAGGTTTCATATGTTTGTGTTCAAAGTAAATATCCTCATTGCTGCAGCTTACTAA
AAATGTAAAGAAAATTGTTTTGCTCGTGTAGATATCTGTAAATTTGTTTTGCATATTAATATATA
TAGATAATTTTTAATAAGCATCCAAGTCTGTTTACTTTAAGAAAACCATTTCCCAAACAGATTTTTTT
TTATTTCAAGAAAATTTGCTACCATTTAAGTAAGAGAAGGTGAGAAGGATGACAGAGGTTGATTGGT
AGCTATTGAATTCATGAAAACCTTTAAGTTAGCATTTGTTAGCAGTTATTATCCAAGCCAGAGTAGGAT
TTGTTACCAGTTGTTATCCAAACCTAATGTTTAAATTACACATTGTTGAAATTAATTACACATTGTTG
ACATGCTTCTCTCCTGATTGTTTTATTTTAACTTGATAGGCATATCTATGAAACCTTTGTAATTT
TAGTTTATTGCTTTACCATTTATTTTACTAGGTAATTAAGAGAACAGATTTTGTCTCTAATTTTTAAG
CCTTATTTACATATGCAGAAACAGCTTAAATATTTTGAAGTATAGACAAACAGTTAATAGATCCACC
ATTAGGAATCAATATATTATGTCATAATAAACATCCTTTTCTTTCACTGAAATTTCTTTAGAAATAA
ACTTATTTTGTGTTGTTATGTTTTGAAACTTGACATAGGATATTTCCCTCTGGCTACACATTCACCTA
CCCTTGTTCTCTATTTAGATTATTCAAATAAAGTTAGTTTGCTTTTATAGTA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_018353.5</u>
Summary:	Required for recruitment of CENPA to centromeres and normal chromosome segregation during mitosis.[UniProtKB/Swiss-Prot Function]
Locus ID:	55320
MW:	37.6