

Product datasheet for SC211231

TSGA10 (NM_182911) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TSGA10 (NM_182911) Human 3' UTR Clone
Symbol:	TSGA10
Synonyms:	CEP4L; CT79; SPGF26
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_182911
Insert Size:	966 bp
Insert Sequence:	>SC211231 3'UTR clone of NM_182911 The sequence shown below is from the reference sequence of NM_182911. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GAAGAGAATCTTTGCTACAGAGATTTCTGACACCTGAAATGATTCTTCACATCCCTGAGAAAGGTCAAA
GTTACAAACTGATTTTTTTTTTTTTGCTACATGAGTGCATTTATCTTTAAATGCTTGGCAATGTTAAA
TGTATTTATTAACCTTGTGTCTCTGAATCTCTGTTCTAATGTCCATGTTGCAGTGATCTGAGATGACT
TATAAAAACAAAATGTATATGGCTCTTTCTATCCATGCAGTGATAGTGAGTGAAAATCTGCTTACTT
CACTATTGAACACTATTCTGTTCACTATCTGGAGTAAATAAAGAAGCTTATTTAAACAGGGAAAAGTGT
TTTTACAAAACCTCTTTCTTTCTTTCTTGGTATCTCAAATAATTGGTTGGTTAATTTTTTTCTTTG
TGATTTGTGTTTTCATGGCTGGAAAATCTGTGCCAAAATGGGTACAACCTACTGATACTGGTCTGGGAT
CTGATATTCCATTTTGAACCAACCTTACATGAATCCATGCCATAATCTGATATTGAAGCAGCAGAAA
GTAGGAGTGTTGTCAGGGTTTTAAACCTATAGTATTGTACCACTCAAAGACTAAGTAACATACATATT
CTACAATTGTGTTTACTACTACTGCTTTTAGAAAGGTAGGTGATGGTATATCTGTCAAGAATGTAAGAT
TTTATTTAAAAATCTTTATTTTAAATAAGTTTTCTTTCAAGTACTGCAATCAACTTAATTGCGATTTTAT
TTTATTGTCTGTGGCAGCTGTCTGTAATCTGGAGTAATAAGATAAGTGATGTTTCCTTCTCCCAAT
TTCTAAATCTCTTACCTTATCTATCCATAGGCTTTGATGTAATGTTTATATTAGAGCCTGAATGCTAC
TGTATTTTACTGATATAAACTATTTGATTCTTGCCATTTTAATAAAAAAAGGCTTGTCATTAATAGTA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_182911.4</u>
Summary:	Plays a role in spermatogenesis (PubMed:28905369). When overexpressed, prevents nuclear localization of HIF1A (By similarity).[UniProtKB/Swiss-Prot Function]
Locus ID:	80705
MW:	37.6