

Product datasheet for **SC211238**

MED15 (NM_015889) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: MED15 (NM_015889) Human 3' UTR Clone
Symbol: MED15
Synonyms: ARC105; CAG7A; CTG7A; PCQAP; TIG-1; TIG1; TNRC7
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_015889
Insert Size: 945 bp
Insert Sequence: >SC211238 3'UTR clone of NM_015889
 The sequence shown below is from the reference sequence of NM_015889. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTCCACCAGGCCCTGCCTCTCAGCCGCCTAGCCAAGACTGCAGGGATGGCCCCGAGCCTCATCGGGGCCA
AGGACACACGCCTCCTGTCAGACACTTCTAGGTGTTGGCTTCTTAGAGAGCCTGGGGTTAGGTTAGCT
TTCCTGCTTTTATCTTCTGCCTTGGGGACCTGCCAAACGAAATCCACACCTGTACAGAAGTGGGATAG
GCGCAGTGGAGCGGGTGTCTGGGGGGCGTTGGCCGACTTCTTAGAGAAGGCCCTCCATGTGACTTCTCCT
CCCAGGAGCCAGATGCGATCCTCAGGCTGCTCTACCGTGGCCTGTCCACGGTCCAGGTCCATCTCAGC
AGCGTGAGGGTGCCTCAGGGTGTGTTAGAGCGTCTCGTGTGTGCTAGACGCACCCCTACTCGTTCCT
ATAGAACACAGAGGACATAGGAAACCCCTAAAACACACATGGGATTCTCTGGTCACAGTTTTGGGTTCA
GGCTATGCTGCTTTGGGCAGGTGGAGCACCCCCGAGGAAGCCTGCAAGTCCAGGGCACAGGCTGCCTT
TTGGAGGGAGGGCTGGCCATAGGTGCTGCTGGCTCCCGCCACCAAGTGGGCTCAGCCCTCACGGCA
TTCCTGCTGAGCACCGTGGGGCACCCAGGGAGCAGGGGCGTCAGGGATCCTGCTGCCGGCACCCCTGTG
CCGCTGGCATGAGGGCCGTGTCCTCCACTGTGAAGGATGAAGAGCAAGGCCCTCAGGACCCGTGCTCTCA
GAGCACACACTGAGCACCCAGAGACAGCGGGCCTGGCAGCGGGCCGGCCATGCAGGGAGCGCCTC
CCTATGTTGCCTGCCACTCTGGGCACCGGCCAGCACCTCTGGTGAGAAGAGGTCCCCCTTTTTATGT
GCACTACCCACCATCTGTGATTATAATAAATTTATTATTCCTGTGTT
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

| | |
|------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | NM_015889.5 |
| Summary: | The protein encoded by this gene is a subunit of the multiprotein complexes PC2 and ARC/DRIP and may function as a transcriptional coactivator in RNA polymerase II transcription. This gene contains stretches of trinucleotide repeats and is located in the chromosome 22 region which is deleted in DiGeorge syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014] |
| Locus ID: | 51586 |
| MW: | 33.7 |