

Product datasheet for SC211247

ZNF439 (NM_152262) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF439 (NM_152262) Human 3' UTR Clone
Symbol:	ZNF439
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152262
Insert Size:	947 bp
Insert Sequence:	>SC211247 3'UTR clone of NM_152262 The sequence shown below is from the reference sequence of NM_152262. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AAGAATGCACTCTGGAGAAAGACCTTAAAATATAAGATATATGGGAAACACTTTTATTCTGCCAAGTT
ATTTCAAACACATGAAAAAATTCACACTGGAGAGAAACCCTATAAATGCAAGCAATGTGGTAAAGCCTT
AATTGTTCCAGTTCCTTTTCGATATCTAAAAGGACTCACAGTGGAGAAAACTCTATGAGTGAAGCAAT
GTGGGAAAGTCCTTCAGATCTGTCAAGAACCTTTCAATTTATGAAAGGACACACACTGGAGAGAAACCCT
ATGAATGTAAGAAATGTGGAAAAGCGTTCCATAATTTCTCTTTTCAAATACATGAAAGTTGCACAG
AGGAGAGGCGCCCTAAGAATGTAAGCATTGTGGGAAAGCATTATCTGCCAAGATCGTTTGAATACA
TGCAAAACACACACTGGAGAGAAACCTATGAATGTAAGGAATGCAACAAGCATTCAATTATTTTCTT
CCTTGATATACATGAAAGGACTCATACGAGAGAGAATCCGTATGAATGTAAGGATTGTGGGAAAGCAT
TCAGCTTGCTTAATTGCTTTCATAGACATGTAAGACACACCAGAAGGAAACCCTATGAATGTAAGCAA
TGTGGCAAAGCTTTCACCTTCTCCAGTCTTTTCAATATCATGAAAGGACTCACACTGGGAGAAACC
GTATCAATGTAAGCAATGTGGGAAAGCCGTCAGATCAGCCTCAAGACTTCAAATGCATGGAAGCACTCA
CACTTGGCAGAACTCTATGAATGTAAGCAGTATGGGAAAGCCTTCAGATCGGCTAGGATTCTTTGAAT
ACAAATGAATGTAACAATTAACCTGTTTATAATACTGTATACTAACAATGTTATTCTTTTAAA
TAATTAAGAAGCTATAATAAAATATCCATTGGTGTGATGATTAGATCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_152262.3
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	90594
MW:	36.9